PDS

Product Data Sheets

IRRIGATION LINE - 1.011.0633

- REFERENCES
- TECHNICAL DATA SHEET
 - Annex 01 – LABELS AND IFU
- COMPONENTS
 - Annex 02 - COMPONENTS
- COMPONENTS FOR PACKAGING
 - Annex 03 – COMPONENTS FOR PACKAGING
- STERILIZATION
 - Annex 04 – STERILIZATION PROCESS VALIDATION
- DECLARATION AND APPLIED STANDARD

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RE-USABLE TUBING FOR CONDUCTING AND PUMPING STERILE LIQUID TO THE IRRIGATION SYSTEM OF A HAND
PIECE OR CONTRA ANGLE.

TECHNICAL DATA SHEET

1. Product identification

Product is classified in Class IIa in accordance with Rule 5 of Annex IX to Directive 93/42/EEC (updated by Directive 2007/47/CE) and are CE marked 0546
GMDN code: 44772

2. General product description and Manufacturer identification

Omnia Spa
Via Francesco Delnevo 190 - 43036 Fidenza – Parma
Tel. 0524.527453 – Fax. 0524.525867 – www.omniaspa.eu
Distributed by:

KaVo Dental GmbH
Bismarcking 39
88400 Biberach/Riß
Tel.: +49 (0) 7351 56-0
Fax: +49 (0) 735156-1488
www.kavo.com

3. Material : See dedicated section "Components"

4. Sterilization: The product is NO STERILE.

Method of sterilization: steam (according to indication in IFU see Annex 01)

Information about the sterilization process validation in Annex 04

5. Packaging

Products are packaged in specific envelope for steam sterilization

	First packaging	Second packaging	Final packaging	Pallet 120x80
No pcs	1	2	22	880
Box /envelope	Plastic/paper envelope Medical paper 60 gr/m ² - Plastic film PET12/PP40 PEEL	White box 25x30x4h cm + N° 1 copy of the IFU	Box 37,5x28x41h cm	/


6. TECHNICAL DRAWING

1.011.0633
(30.F0341)

Components - Technical description:

(see section "Components" for all technical data)

1	2 ways vented spike - PC and ABS
2	Clamp - PP
3	Connector - POM
4	Plastic Key - POM
5	Connector for pump tube - APEC1745 (Bayer) transparent
6	Stainless steel ring
7	Fixing clips - ULTRASON (BASF)
20	Irrigation tube Silicone 60 ShA Ø mm 3x4,1 - 165 mm
21	Tube for peristaltic pump Silicone 50 ShA - Ø mm 4x7 - 190 mm
22	Irrigation tube Silicone 60 ShA Ø mm 1x3 - 2420 mm
23	Tube Silicone 60 ShA Ø mm 4x7 - 30 mm

OMNIA Spa—Via F. Delnevo 190—43036 Fidenza (PR)

pagina 3

7. LABELING

In Annex 01 :

LABEL FOR FIRST PACKAGING - 1 PC
Label real dimensions 10 x 15 cm

LABEL FOR SECONDARY PACKAGING - 2 PCS
Label real dimensions 10 x 15 cm

LABEL FOR FINAL PACKAGING - 22 PC
Label real dimensions 10 x 15 cm

8. INSTRUCTION FOR USE :

In Annex 01 : IFU with all specification and information according to ISO 17664

COMPONENTS

In Annex 02 :

all details of components in point 6 - drawings and technical data sheets of the materials

COMPONENTS FOR PACKAGING

In Annex 03 :

all details of envelope - drawing and technical data sheet

*Envelope for steam sterilization
(Medical paper 60 gr/m² - Plastic film PET12/PP40)
20X30 cm*

10.06.2015 - Annex 1

ANNEX 1

LABELS and IFU

First packaging	Plastic/paper envelope 20 x 30 cm Medical paper 60 gr/m ² - Plastic film PET12/PP40 PEEL Label 10 x 15 cm
------------------------	--

Barcode construction

	LIC	Ref.	Packaging level		Q.ty field		Lot field		Check digit	
BLISTER (1 PC)	+EFID	1.011.0633	0	/\$\$	8		7	<assigned> e.g. 123456	<calculated> 1	
BOX (2 PCS)			1						02	<calculated> 3
CARTON (22 PCS)			2						22	<calculated> 6

Distributed by:
KaVo Dental GmbH
Bismarcking 39
88400 Biberach/Riß
Tel.: +49 (0) 7351 56-0
Fax: +49 (0) 735156-1488
www.kavo.com

Sterilisierbares Schlauchset S600
Sterilizable Hose Set S600

Catalog number NOT MADE WITH NATURAL RUBBER LATEX
REF 1.011.0633 PHTHALATES FREE

Lot number
LOT

Date of manufacture

 Manufacturer: OMNIA Spa
www.omniaspa.eu **CE** 0546

Consult instruction for use

Keep away from sunlight

Keep dry

Do not use if package is damaged


Non-sterile

! Caution: US Federal law restricts this device to sale by or on the order of a health care professional/dentist. For dental use only.

Rev. 01/2015

Second packaging	Carton box - white - 25x30x4h cm - 2 pc + IFU Label 10 x 15 cm
-------------------------	---

Distributed by:
KaVo Dental GmbH
Bismarcking 39
88400 Biberach/Riß
Tel.: +49 (0) 7351 56-0
Fax: +49 (0) 7351 56-1488
www.kavo.com

Sterilisierbares Schlauchset S600
Sterilizable Hose Set S600
BOX =2 PCS

NOT MADE WITH NATURAL RUBBER LATEX
PHTHALATES FREE

Catalog number
REF 1.011.0633

Lot number
LOT

Date of manufacture

Manufacturer:
OMNIA Spa
www.omniaspa.eu

CE 0546

Consult instruction for use 

Keep away from sunlight 

Keep dry 

Do not use if package is damaged 

Non-sterile 

Caution: US Federal law restricts this device to sale by or on the order of a health care professional/dentist. For dental use only.

Rev. 01/2015

Final packaging	Carton box – white – 37,5x28x41h cm – 22 pc Label 10 x 15 cm
------------------------	---

Distributed by:
KaVo Dental GmbH
Bismarckring 39
88400 Biberach/Riß
Tel.: +49 (0) 7351 56-0
Fax: +49 (0) 735156-1488
www.kavo.com

Sterilisierbares Schlauchset S600

Sterilizable Hose Set S600

11 BOX =22 PCS

NOT MADE WITH NATURAL RUBBER LATEX
PHTHALATES FREE

Catalog number
REF 1.011.0633

Lot number
LOT

Date of manufacture

Manufacturer:
OMNIA Spa
www.omniaspa.eu

CE 0546

Consult instruction for use

Keep away from sunlight

Keep dry

Do not use if package is damaged

Non-sterile


! Caution: US Federal law restricts this device to sale by or on the order of a health care professional/dentist. For dental use only.

Rev. 01/2015

IFU – 80.M0050 – REV00/APRILE/2015

- **Paper 60 gr A5 format**
- **7 languages**
 - **Italian**
 - **English**
 - **German**
 - **Russian**
 - **Japan**
 - **Chinese**
 - **Spanish**

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ANNEX 2

COMPONENTS

COMPONENT 1 : 2 WAYS VENTED SPIKE—PC and ABS

PHYSICAL FEATURES	NORMS	TYPICAL VALUE	MEASURE UNITS
1) Specific weight	ISO 1183	1.2	Kg/dm ³
2) Hardness	ASTM D-2240	D/81±3	SHORE
3) Tensile strength	ISO 527	63 at yield	MPa
4) Elongation	ISO 527	6 at yield	%
5) Thermal coefficient of linear expansion	ASTM E-831	0.7	10 ⁻⁴ /K
6) Stiffness in torsion			
7) Vicat softening point			
STABILIZER(S)	:		
PLASTICIZER(S)	:		
COMPLIANCE	:	-Meets the requirements of: U.S.P. XXIII Ed. - Class VI - Plastics	
MAJOR USE	:	Extrusion <input type="checkbox"/>	
	:	Moulding <input checked="" type="checkbox"/>	
	:	Blowing <input type="checkbox"/>	
PACKAGING	:	Goods delivered in polyethylene or paper bags weighing 25 kg each or 600-800 kg	
ASPECT	:	Plastic granules	Colourless <input checked="" type="checkbox"/>
			Red <input type="checkbox"/>
			Blue <input type="checkbox"/>
			White <input type="checkbox"/>
STORAGE CONDITIONS	:	Maintain containers closed in a clean and well-aired environment. Avoid direct light and excessive humidity.	
STERILIZATION TYPE	:	Ethylene oxide <input checked="" type="checkbox"/>	
		Gamma radiation <input type="checkbox"/>	
		Steam <input checked="" type="checkbox"/>	
NOTES	:	-This material can be used for components which come into contact with blood and/or "solutions".	

COMPONENT 2: CLAMP —PP

Product Data and Technical Information

[✉ Email this page](#) [🖨 Print this page](#)



Name Moplen

Grade HP501H

Resin type Polypropylene, Homopolymer

Description *Moplen* HP501H is a homopolymer for injection moulding which exhibits a good stiffness/impact balance at ambient temperatures. The main applications of *Moplen* HP501H are caps, closures and furniture. *Moplen* HP501H is suitable for food contact.

COMPONENT 3 : CONNECTOR FOR PUMP TUBE —POM


KOCETAL[®] K300

Kocetal K300 is a general injection molding grade with a good mechanical property for wide range.

Property	Test Method	Units	Value
Physical			
Melt Index	ASTM D1238	g/10min	9.0
Specific Gravity	ASTM D792	-	1.41
Shrinkage	ASTM D955	%	2.0
Water Absorption 23°C Equilibrium 60%RH	ASTM D570	%	0.22
Thermal			
Melting Point	ASTM D1525	°C	166
Heat Deflection Temperature under Load			
HDT/A 0.45 MPa (4.6 kg/cm ²)	ASTM D648	°C	158
HDT/A 1.82 MPa (18.6 kg/cm ²)		°C	110
Vicat Softening Temperature	ASTM D1525	°C	162
Coeff. Of Linear Thermal Expansion	ASTM D696	× 10 ⁻⁵ cm/cm·°C	13
Flammability	UL 94		HB
Mechanical			
Tensile Strength 23°C	ASTM D638	MPa (kg/cm ²)	64 (650)
Tensile Elongation 23°C	ASTM D638	%	60
Flexural Strength 23°C	ASTM D790	MPa (kg/cm ²)	94 (960)
Flexural Modulus 23°C	ASTM D790	MPa (kg/cm ²)	2,500 (25,500)
Notched Izod Impact Strength 23°C	ASTM D256	J/m (kg-cm/cm)	67 (6.8)
Rockwell Hardness M scale	ASTM D785	-	80
Electrical			
Surface Resistivity	ASTM D257	× 10 ¹⁶ Ω	1
Volume Resistivity	ASTM D257	× 10 ¹⁴ Ω-cm	1
Dielectric Strength	ASTM D149	kV/mm	19
Dielectric Constant	ASTM D150	-	3.7
Dielectric Dissipation FactorConstan 10 ⁶ Hz	ASTM D150	-	0.006

COMPONENT 4 : PLASTIC KEY —POM - Color: white RAL 9016

Based on KaVo Dwg No 1.009.8819

Certification

Company : **KOLONPLASTICS, Inc.**
 Address : **1018 Ungmyung-Dong, Gimchon-Si, GyungSangBuk-Do, 740-180, Korea**
 Name and the function :

We certify the co-polymers, **POM KOCETAL® K300; 1,3,5-Trioxane, polymer with 1,3-dioxolane, Cas No.24969-26-4 (submitter names) : copolymer of dioxolane and trioxane** , supplied by **KOLONPLASTICS, Inc. does not contain Environmentally Hazardous Substances mentioned category.**

- **Materials for Certification**
POM KOCETAL® K300

- **Category**
**PARABENS
 PHTHALATE
 alkylphenols**

13 July 2011
 Date

R&D Division
 KOLONPLASTICS, Inc.
 Tel: +82-54-420-8472
 Fax: +82-54-420-8329

Approved by 
Jong-Mun Kim / R&D Division Manager
 Signature & Supplier's stamp

COMPONENT 5 : CONNECTOR

Apec 1745

Standard grades / Medical applications, suitable for superheated steam sterilization


Easy to demold, suitable for superheated steam sterilisation up to 143 °C as well as for pharmaceutical applications according to United States Pharmacopeia (USP) XXII Class VI, softening temperature (VST/B 120)=170 °C

COMPONENT 6 : STAINLESS STEEL RING

COMPONENT 7 : FIXING CLIPS —ULTRASON (BASF) – PSU- TRANSPARENT

Based on KaVo Dwg No 1.009.8818

Product Information

July 2004

**Ultrason[®] S 2010 natural
(PSU)**

BASF
The Chemical Company

Product description

Unreinforced, medium viscosity standard injection moulding grade.
Abbreviated designation according to ISO 1043-1: PSU

COMPONENT 20 : SILICONE 60 ShA - 3 x 4,1—L 16,5 cm
COMPONENT 22 : SILICONE 60 ShA - 1 x 3 — L 242 cm
COMPONENT 23 : SILICONE 60 ShA - 4 x 7 —L 3 cm

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ANNEX 3

COMPONENTS of PACKAGING

Envelope for steam sterilization (Medical paper 60 gr/m² - Plastic film PET12/PP40) – 20 x 30 cm

Property chemical-physical	Standards	Value	Unit of measure
PET 12μ			
Density	ASTM D1505	1,395 + 1,405	g/cm ³
COF static	ASTM D1894	< 0,6	/
COF kinetic	ASTM D1894	< 0,55	/
Point/fusion interval	/	260	°C
PP 40μ			
Density	Internal Method	0,9	g/m ³
COF	ASTM D1894	0,15 + 0,30	/
Point/fusion interval	/	150 + 170	°C
ADHESIVE			
Weight	/	1,6 – 2	g/m ²
<p>Test made in environmental conditions conform to Standard UNI EN ISO 11607-1:2006 (4.4.3). (23 ± 1) °C – (50 ± 2) % UR – 24 h, a part from different conditions indicated on the applied Standards.</p>			
CONFORMITY TO STANDARD	<p>Composite Film PET12/PP40 green is conform to Standard UNI EN ISO 11607:2006 and to Directive 94/62/CE of European Parliament and Council, , dtd December 20, 1994; on packaging and packaging' s waste.</p> <p>Conformity of the components to Standards applied:</p> <p>PET POLYETYLENE – UE: Commission Directive 2002/72/CE, 2004/19/CE, 2005/79/CE, 94/62/CE, 2002/16/CE, 2007/19/CE, Regulation 683/98, 1935/2004 and 1895/2005; Italy: Ministerial Decrees 21.3.1973 and DPR 777/82; USA: FDA, CFR, title 21 § 177.1500 (a) (6) e (b) (6.1-6.2) e CFR, title 21 § 177.1630</p> <p>PP POLYPROPYLENE – UE: Commission Directive 2002/72/CE; Italy: Ministerial Decrees 26.4.1993 n. 220, 28.10.1994 n. 735, 24.9.1996 n. 572, 6.2.1997 n. 91, 22.7.1998 n. 338, 15.6.2000 n. 210, and Ministerial Decrees 21.3.1973 n. 104, attach II, section I, part B; USA: FDA, CFR, title 21 (2002) § 177.1520 (a) (1), (b) and (c) 1.1</p> <p>ADHESIVE – UE: Commission Directive 2002/72/CE; Italy: Ministerial Decrees 21.3.1973; USA – FDA, CFR, title 21 § 175.105</p>		
DESCRIPTION MATERIAL	<p>Composite Film PET12/PP40 green is composed of a laminated Polyester and a laminate of Polypropylene. The adhesive interposes of two Films does not contain solvent, so the final product is solvent free.</p>		
APPLICATION	<p>Composite Film PET12/PP40 green is especially developed for single use medical application ; suggested for rotating or flat welding. Composite Film could be weld on itself and paper support. This material is suitable for <i>Steam</i> and EO (gas) sterilization.</p>		
TRACEABILITY	<p>The product/raw material is guaranteed through documentations keep by QC Encaplast's files for a period of 10 years.</p>		
GENERAL CONDITIONS OF STORAGE	<p>15-25°C 60% UR</p>		
 <p><i>It's responsibility of the user to validate the single packaging in accordance with Enclosure I Medical Directives Device 93/42/CEE, since product supplied by Encaplast is considered a component of Medical Device.</i></p>			
Property of ENCAPLAST Re-production strictly subject to authorization.			Drawer:  Date: 31-03-08 Variation N°: 256/08
	Responsible:  Date: 31-03-08 Date: 31-03-08		
	Description Composite Film PET 12 / PP 40 Green (Data Sheet)		
	Code: 2503350000	Rev. C	Pag. ___ of ___

OMN

PROPERTY PHYSICAL-CHEMICAL	Standards	Typical Value	Unit of Measure
Paper Grammage	ISO 536	60	g/m ²
pH water extract	ISO 6588-2	7	/
Chloride content	ISO 9197	< 0.05	%
Sulphate content	ISO 9198	< 0.25	%
Burst strenght (dry)	ISO 2758	350	KPa
Burst strenght (wet)	ISO 3689	150	KPa
Traction strenght long. (dry)	ISO 1924-2	6,4	kN/m
Traction strenght transv.(dry)	ISO 1924-2	3,4	kN/m
Traction strenght long. (wet)	ISO 3781	2,1	kN/m
Traction strenght transv. (wet)	ISO 3781	1,1	kN/m
Tear strenght long.	ISO 1974	600	mN
Tear strenght transv.	ISO 1974	650	mN
Air permeability	ISO 5636-3	11,4	µm/(Pa.s)
Superficial roughness Bendtsen FS	ISO 8791-2	375	ml/mn
Superficial roughness Bendtsen WS	ISO 8791-2	300	ml/mn
Superficial absorption (Cobb)	EN 20535	15	g/m ²
Hydrophobia (water repellency)	ASTM D779-03	35	sec
Fluorescence	DfN 58953-6	0	%
Pore diameter	EN 868-2 (all. C)	21	µm

Test made in environmental conditions conform to standard UNI EN ISO 11607-1:2006 (4.4.3), (23 ± 1) °C – (50 ± 2) % UR – 24 h, unless different conditions indicated on specific standards applicated.

CONFORMITY TO STANDARDS Medical Paper HL 60g/m² is conform to standards :
UNI EN ISO 11607-1:2006, EN 868-3:1999, FDA CFR 21. BfR (B_gVV) XXXVI;
and is made in accordance to standard EN ISO 9001.

MATERIAL DESCRIPTION Medical Paper HL 60g/m² is humidity-resistance and have a degree of roughness that made it suitable to EO (gas) and Steam sterilization. It present high bacterial barrier. Is indicated as a support for flex graphic and rotogravure printing.

APPLICATION Medical Paper HL 60g/m² is especially indicated for the production of pouches and rolls, printing and unprinting, where the Paper is hot welding with composít film polyester/poliprophylene. This packaging can be used for EO (gas) and Steam sterilization. If is welding with polyester/polyethylene, is indicated for EO (gas) and Radiation.

TRACEABILITY Traceability of products/prime material is guaranteed throught registration documents held in CQ Encaplast file, for a period of 10 years.

GENERAL CONDITION OF STORAGE 15-25°C 60% UR

It's the responsibility to the user to validate the single packaging in accordance with the Enclosure 1 Medical Device Directives 93/42/CEE, since product supplied by Encaplast is considered a component of the Medical Device.

Property of ENCAPLAST Re-production strictly subject to authorization.		Drawer <i>Marchi A.</i>	Date 31-03-08	Variation N° 256/08
		Responsable <i>Neri Maria</i>	Date 31-03-08	Date 31-03-08
Description MEDICAL PAPER HL 60g (Material Sheet)				
Code 2501000000		Rev. D	Pag. <u> 1 </u> of <u> 1 </u>	

OMN

10.06.2015 - Annex 4

ANNEX 4

STERILIZATION PROCESS VALIDATION

The sterilizable KaVo hose set is comparable, for components, intended use and mode of use to the sterilizable line used for sterilization process validation, which defines:

1. Sterilisation of the line at 121 ° C by steam (autoclave tipoB)
2. Re-sterilizable for a maximum of 25 times
3. Method of reconditioning of the line after each use

The reports show that the line does not alter the components and functionality as a result of the 25 cycles, as it is not altered the functionality.

These are the standards applicable for protocol validation:

- European Pharmacopoeia – current ed.
- UNI EN ISO 11737-1:2006 → Sterilization of medical devices – Microbiological methods – Part 1: Determination of population of microorganisms on products
- UNI EN ISO 11737-2:2010 → Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
- UNI EN ISO 11138-1:2006 → Sterilization of health care products – Biological indicators – Part 1: General
- UNI EN ISO 11138-3:2009 → Sterilization of health care products – Biological indicators – Part 3: Biological indicators for moist heat sterilization process

Report 14/2011 – 08/08/2011 → Validation of steam sterilization

Report 39/2012 – 27/07/2012 → Technical assessment

Report 40/2012 – 27/07/2012 → Evaluation of toxicity

For each production, we carry out the bioburden analysis in order to define and monitor the level of contamination, in accordance with the limits identified in the validation of the sterilization process.

Specifically, on the KaVo hose set, tests were conducted to evaluate the re-sterilization of the seal assembly of the component 1 with component 20 using the silicone tube component 23. Functionality tests were performed and confirmed by the customer on the samples provided with lot 15150.