



An evidence-based technique
option for non-aerosol
generating clinical practice^{1,2}



In today's environment, there is a focus on restorative techniques where aerosols are reduced to a minimum. One important non-aerosol generating technique option is called **Atraumatic Restorative Treatment (ART)**.

ART is a minimally invasive operative intervention where hand instruments alone are used to prepare a cavity atraumatically. The cavity is subsequently restored with a self-adhesive, bulk fill, high strength glass ionomer restorative. The basis for this technique comes from numerous clinical studies on partial removal of carious tissue followed by sealing and restoration with glass ionomer cement.¹⁻⁴

According to research evidence and published recommendations on caries removal, a significant outcome is the recognition that only soft, decomposed tooth tissue needs to be removed and the remaining softened, partly demineralised dentine on the pulpal floor can be left behind to remineralise.¹⁻⁴

As well as the benefits for the tooth of a minimal intervention approach, the ART technique can offer dental professionals a simplified, yet evidence based system^{1,2} that is less traumatic for the patient, making it an ideal option for today's challenging climate.

Overview

The origin of ART technique comes from research in the 1990s to develop a restorative and preventive sealant technique suitable for developing economies where access to electricity and cavity cutting equipment was limited and an alternative to extraction was needed.^{1,2}

The original research was a collaboration with the World Health Organisation and GC became active participants after developing a new high strength glass ionomer cement for ART called Fuji® IX. Fuji® IX was so called due to the 9-steps taught by WHO as basis for the ART technique.

Products

ART instruments for cavity preparation

- Enamel access cutter
- Hatchet (to improve access to the cavity if required)
- Spoon excavators
- High strength posterior glass ionomer cement, e.g. GC Fuji® IX GP EXTRA, GC Fuji® BULK.
- In more aesthetic, non-stress bearing areas a resin reinforced light cured glass ionomer restorative could be used e.g. GC Fuji® II LC
- A dentine conditioner is recommended prior to placement
- A transitional coating is recommended if moisture control during initial set is difficult e.g. Vaseline or cocoa butter.

More recently a modification of this technique has been promoted called Simplified Modified ART (SMART) where access to electricity has meant a capsulated glass ionomer is used and newer, even less invasive cavity preparation principles have been adopted.

An alternative SMART (Silver modified ART) technique is actively promoted in USA where the use of Silver Diamine Fluoride is advocated as a pre-treatment prior to glass ionomer cement application to restore a cavity.



GC Fuji® IX GP EXTRA



GC Fuji® BULK



GC Fuji® VII EP*



DENTIN CONDITIONER



COCOA BUTTER

Technique

- Isolation
- Remove plaque with a wet cotton wool pellet, then dry with a dry pellet
- If required, widen the entrance to the cavity using an enamel access cutter or dental hatchet. Place the blade or tip of the instrument into the cavity and turn forward and backward like turning a key in a lock. This movement chips off small pieces of carious enamel.
- Carious dentine is removed with the excavators using circular scooping movements. It is important to remove all the soft caries from the enamel-dentine junction
- Clean the cavity with wet and then dry cotton wool pellets.
- Apply CAVITY CONDITIONER for 10 seconds to both dentine and enamel. Wash the cavity with a wet cotton wool pellet at least twice, then dry with a dry pellet
- Activate and mix Fuji BULKcapsule (or alternative) for 10 seconds
- Immediately inject the mixed cement direct into the cavity and adapt
- If using a “press-finger technique”, rub a small amount of petroleum jelly on the gloved index finger and cement firmly into the cavity and fissures. Remove the finger sideways after a few seconds. The time from the start of mixing until removal of finger should be no more than one minute. Remove excess
- Check the bite. Use the blade of a carver to make occlusal adjustments



1. Place cotton roll at buccal area.



2. Remove soft carious lesion at DEJ by spoon excavator in circular motion only (partial caries removal).



3. Clean floor of the cavity with wet and dry cotton pellets.



4. Apply DENTIN CONDITIONER using a cotton pellet for 20 seconds to remove the smear layer. Clean with wet and dry cotton pellets afterward.



5. Press finger on top of the filling then slide finger out carefully. Carve the excess of GI with spoon excavator. Avoid moisture contamination and dry-out.



6. Coat GI with COCOA BUTTER or Vaseline. Check occlusion by biting from opposite plane and carve the excess GIC.

ART Technique for 2020 and Beyond

An atraumatic treatment strategy for **anxious patients**, especially useful for initial carious lesions

Core to oral health strategies for the **elderly and dependent** in home care and residential care facilities

A "best practice" strategy for treating **deep carious** lesions

A treatment strategy for **non-aerosol** generating restorative procedures

A treatment option that can **reduce** the need for restorative care under **GA and sedation**

ART sealants targeting **erupting molars** are an important public health strategy

Evidence

ART technique:

1.de Amorim RG, Frencken JE, Raggio DP, Chen X, Hu X, Leal SC. Survival percentages of atraumatic restorative treatment (ART) restorations and sealants in posterior teeth: an updated systematic review and meta-analysis. Clin Oral Investig. 2018; 22:2703-2725.

2.Frencken JE. Atraumatic restorative treatment and minimal intervention dentistry. Br Dent J. 2017; 223:183-189.

Partial caries removal:

3.Bannerjee A, Frencken JE, Schwendicke F, Innes NPT. Contemporary operative caries management: consensus recommendations on minimally invasive caries removal. Br Dent J 2017; 223:215-222

4.Ngo H C, Mount G, Mc Intyre J, Tuisuva J, Von Doussa RJ. Chemical exchange between glass-ionomer restorations and residual carious dentine in permanent molars: An in vivo study. J Dent 2006; 34:608-613.

DISCLAIMERS:

The products cited in this leaflet are for professional use only. Always read and follow the instructions for use.

Descriptions: DENTIN CONDITIONER is a mild polyacrylic acid solution designed to remove the dentinal smear layer and to condition dentine and enamel before the application of your glass ionomer restorative. CAVITY CONDITIONER is a 20% polyacrylic acid solution designed to remove the dentinal smear layer and to condition dentine and enamel before the application of your glass ionomer restorative. It will increase the bond between your glass ionomer cement and tooth structure for added longevity. GC Fuji IX GP EXTRA and GC Fuji BULK are multipurpose auto-cure glass ionomer cements used for restoring teeth. GC Fuji VII EP is a glass ionomer cement with CPP-ACP for tooth protection. Once applied, GC Fuji VII EP forms a mechanical barrier to protect the covered and surrounding tooth surfaces and to minimise plaque accumulation. COCOA BUTTER is a separating medium and protective coating for the prevention of water and saliva contamination to exposed surfaces of glass ionomer cements during the first 24 hours after placement. GC Fuji II LC is light-cured glass ionomer restorative used for restoring teeth. GC Fuji IX is a glass ionomer restorative used for restoring teeth.

Contraindications: The products cited are contraindicated for pulp capping. Avoid use of the cited glass ionomer range of products in patients with known allergies to glass ionomer cement. Avoid the use of DENTIN CONDITIONER, CAVITY CONDITIONER in patients with known allergies to polyacrylic acid. Do not use GC Fuji VII EP in patients with a proven or suspected milk protein allergy or sensitivity and in patients with known allergies to glass ionomer cement. Avoid use of GC Fuji II LC in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer. Avoid the use of COCOA BUTTER in patients with known allergies to cocoa butter or nut. Avoid use of Fuji BULK in patients with known allergies to glass ionomer cement. Avoid use of Fuji IX in patients with known allergies to glass ionomer cement. In rare cases, the products cited in this leaflet may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the products and refer to a physician.

Cautions: DENTIN CONDITIONER and CAVITY CONDITIONER: In case of contact with oral tissue or skin, flush immediately with water. In case of contact with eyes, flush immediately with water and seek medical attention. GC Fuji BULK, GC Fuji VII EP, GC Fuji II LC, GC Fuji IX GP EXTRA: In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. In case of contact with eyes, flush immediately with water and seek medical attention. COCOA BUTTER: In case of contact with eyes, flush immediately with water and seek medical attention.



Since 1921
Towards Century of Health

GC Australasia Dental Pty Ltd
1753 Botany Road
Banksmeadow NSW 2019
Australia

T: +61 2 9301 8200
E: info.australasia@gc.dental
www.gcaustralasia.com