B.A. INTERNATIONAL

Ultimate BASE280

LED Curing Light

INSTRUCTIONS FOR USE

Lámpara LED de fotocurado

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MODEL: BASE280

P/N: BA110150

Lampe à Polymériser LED

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LED-Polymerisationslampe

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Lampada polimerizzante a LED

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Led-uithardingslampmp

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MADE IN TAIWAN
Fabricado en Taiwan
Fabriqué en Taiwan
Hergestellt in Taiwan
Prodotto in Taiwan
Gemaakt in Taiwan

B.A. INTERNATIONAL LTD.
Unit 9, Kingsthorpe Business
Centre Studland Road,
Northampton - NN2 6NE,
U.K.
Rev. 2019/04

Disclaimer:

1. First of all, we would like to thank you for purchasing BA Ultimate BASE280 LED Curing Light. Please carefully read the user manual and operate the device accordingly. Please keep this user manual well preserved for future reference.
2. BA International shall not assume any responsibility if any malfunction, damage, or accident resulting in bodily injury caused by improper removal, modification, maintenance, or repair of the product conducted by unauthorized personnel.
3. BA International shall not assume any responsibility if any malfunction, damage, or accident resulting in bodily injury is caused by the use of this product in conjunction with other regulated products.
4. BA International shall not assume any responsibility if any malfunction, damage, or accident resulting in bodily injury is caused by improper removal, modification, maintenance, or repair of the product with components not supplied by the manufacturer.
5. BA International shall not assume any responsibility if any malfunction, damage, or accident of bodily harm is caused by improper operation that is incompliant with the manual instructions.
6. BA International shall not assume any responsibility if any malfunction, damage, or accident resulting in bodily harm is caused by inappropriate modifications to the power supply, installation environment, or any other safety-regulated operating conditions as specified in this manual.
7. BA International shall not assume any responsibility if any malfunction, damage, or accident resulting in bodily harm is caused by natural disasters, environmental causes, or acts of God.

Content:

• Intended Use.....	3
• Safety precautions.....	4
• Explanation of Graphical Symbols.....	3
• Parts & Functions.....	4
• Installation.....	4
• Panel operation.....	4
• Charger signals.....	8
• Troubleshooting.....	9
• Specifications.....	9
• EMC Tables.....	9

Intended Use:

This “Ultimate BASE280 LED Curing Light” is a dental curing unit utilizing visible light programmed for the polymerization of light curing materials in the wavelength range of 360-420nm (UV Light) and 420-490nm (Blue Light) by dental professionals.

Safety Precautions:

Please follow the instructions in this manual for correct and safe operation of the device. Pay special attention to the following warning signs beside the operation descriptions where applicable.

DANGER :

This message appears where the operation may potentially cause severe injury if instructions are not followed correctly.

WARNING :

This message appears where improper operation may cause serious damage or defect to other objects if instructions are not followed correctly.

CAUTION :

This message appears where improper operation may cause slight damage to other objects or the human body if instructions are not followed correctly.

IMPORTANT :

This message appears where the operation may involve risk of defecting the device itself if not properly conducted.

EU-wide legislation, as implemented in each Member State, requires that waste electrical and electronic products carrying the mark (left) must be disposed of separately from normal household waste. This includes monitors and electrical accessories, such as signal cables or power cords. When you need to dispose of your display products, please follow the guidance of your local authority, or ask the shop where you purchased the product, or if applicable, follow any agreements made between yourself.

The mark on electrical and electronic products only applies to the current European Union Member States.

IMPORTANT :

If light curing materials not within the range of 360-420nm (UV Light) and 420-490nm (Blue Light) is used with the device, the material may have difficulty completely curing or required an extended curing duration.

Explanation of Graphical Symbols :

IEC 60601-1 : Applied part classification - Complies with standards protecting against electric shock for type B equipment

Class II equipment

ISO 7000-1641 : Follow operating instructions or consult instructions for use

IEC 60417-5009 : STAND-BY

IEC 60417-5031 : Direct current.

Manufacture date

Manufacturer

Storage conditions: -20°C to 65°C

Cooling System :

- >> This device is designed with a built-in cooling fan which automatically activates whenever the device is in use, and continues to run for a few seconds after the cycle has ended.
- >> If the temperature gets exceeds 45°C, the fan will continue to run until the internal temperature has cooled down.
- >> If the temperature gets exceeds 65°C, the device will automatically shut off and the fan will continue to run. Until the internal temperature has cooled down, the device will be unable to start.
- >> If the unit needs to be used for continuous curing, please decrease the curing time.
- >> Under the premise of not touching the oral tissue, when using the curing light, switching the treatment area or switching any operation settings, it is necessary to cool down for 30 seconds before starting the next curing cycle.

Sterilization of Fiber-Optic Probe :

After each treatment, the fiber-optic probe should be disinfected or sterilized. Wrap completely in an autoclave sterilization pouch, then autoclave at 132°C or 134°C/75 psi for 4 minutes by using a pre-vacuum sterilizer to get maximum results. Allow the fiber-optic probe to air-dry for 20-30 minutes in the sterilization pouch before removing.

IMPORTANT :

Disconnect the fiber-optic probe from the handpiece and autoclave by itself.

IMPORTANT :

Every component of this device is reusable, but only the fiber-optic probe needs to be disinfected or sterilized. Other components can not be disinfected or sterilized.

Product Cleaning:

If there is any dirt or stains on the surface of control panel, handpiece or other components, please use a cloth to wipe clean. If it is hard to clean, please use a neutral cleansing agent. 70% isopropyl alcohol is recommended.


IMPORTANT :

Never use a solvent or volatile oil for cleaning; this can lead to fading of the print on the casing or the metal ring falling off from the probe.


Parts & Functions:

PART A The LED Main Unit :

- (1) Fiber-Optic Probe
- (2) Display Panel
- (3) Mode Button (▲)
- (4) Time Button (▼)
- (5) Li-Ion Battery
- (6) Charging Port
- (7) ⏪/Set Button
- (8) Charging Indicator
- (9) Ready Indicator
- (10) Built-in Sensor
- (11) Display Panel
- (12) °C/°F Switch

PART B Accessories :

WARNING :

Do not modify this device without authorization of the manufacturer.

IMPORTANT :

Do not use other adapters without authorization from the manufacturer.

CAUTION :

Please check the specification before using the adapter, the incorrectness could cause the malfunction of the product.

CAUTION :

Please use only the Ø10 mm fiber optic probe provided by the manufacturer.

Installation:

Please carefully follow the instructions below to ensure proper function of the device.

Fiber-Optic Probe Installation / Replacement

1. Insert the fiber-optic probe into the handpiece until the base of the probe is completely matched to the handpiece and listen for a “click” that signifies the probe has been securely fastened.
2. Place the protection shield on the front end of the handpiece, and the orange anti-glare cone or tube onto the tip of the fiber-optic probe and adjust its position according to the dentist's preference.
3. To remove, simply pull the fiber-optic probe outwards to remove it from the handpiece.

IMPORTANT :

Don't use excessive force when installing the above parts to avoid possible breakage or damage.

Battery Installation / Replacement

IMPORTANT :

After battery replacement, all user settings will revert back to the factory settings.

IMPORTANT :

Prior to using a new battery, please charge for a minimum of 4 hours to ensure maximum battery life.

WARNING :

There is a danger of explosion if the battery is incorrectly replaced. Trained service personnel should replace only with the same or equivalent battery type as recommended by the manufacturer.

Panel Operations:

(A) Power On

Press the “ \odot /Set” button to turn on the display panel. Press the “ \odot /Set” button again to begin curing.

Press again to stop the device from emitting in the middle of a curing cycle.

Note: The unit will automatically snooze after 30 seconds of idle time to preserve battery power.

Note: The use of this device does not require any direct contact with patients.

Note: It is recommended to use disposable sleeves before operating the device.

Note: This device is only to be used by trained personnel in dental practices.

WARNING :

When the device is in use, do not look directly into the light beam. Looking directly into the light beam for an extended duration will cause eye damage.

CAUTION :

Avoid direct contact of light guide with skin, oral tissue and tooth surface. Before use, ensure to put anti-glare cone on the tip of the light guide and keep the light guide at a distance of minimum 2mm from oral tissue or tooth surface.

CAUTION :

It is forbidden for the Fiber Optic Light Guide to be directly in contact with any of the teeth surface or any other oral tissue.

CAUTION :

After each curing cycle, it is necessary to cool down for more than 30 seconds before starting the next cycle.

Any clinical diagnosis need to be enforced by doctor according to the curing environment and each patients' situation.

(B) Select the Mode

Press the “Mode(\blacktriangle)” button to select the mode desired.

M1: Full Power Mode - Full power light intensity.

M2: Ramp-Up Mode - Light intensity gradually climbs from low to high.

M3: Pulse Mode - Alternates between strong and weak light intensity every half second.

Note: When replacing or reinstalling the battery, the default operation mode of the device is “M1”.

Note: When the device snoozes, it will restart with the last mode selected.

(C) Select the Curing Time

Press the “Time(\blacktriangledown)” button to select the curing time desired.

M1: Full Power Mode(1, 2, 3 secs)

M2: Ramp-Up Mode(5, 10 secs)

M3: Pulse Mode (10, 20, 30 secs)

Note: Default curing cycle for M1 is 3 seconds, M2 is 5 seconds, and M3 is 10 seconds.

Note: After replacing the battery, all user operation settings will revert back to the factory settings.

Note: To avoid overheating, it is recommended to avoid overextending the cycle duration.

WARNING :

Heat development (burn hazard)

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Pro-longed exposure of areas near the pulp and soft tissues may result in irreversible or reversible damage.

Therefore, this high-performance curing light must only be operated by trained professionals.

CAUTION :

When the curing light being used for over 10 seconds, it is recommended to cool down for 5 minutes before starting the next curing cycle to avoid the patient being burned by the high temperature fiber-optic probe.

(D) Using the Built-in Radiometer

1. Ensure the charger stand is plugged in.
2. Press the fiber-optic probe flat onto the radiometer sensor.
3. Start a curing cycle and take the light intensity reading.

Note: Any movement of the device while measuring the light intensity will result in variations of the reading.


Note: To Change the temperature reading measurement unit, press the °C/°F switch.

(E) Handpiece Display Panel Descriptions

- M1** : Full Power Mode, the most using in clinical practice - to boost with the full power light output for ultrafast performance to cure all resin-based restoratives in 1, 2 or 3 seconds.
- M2** : Ramp-Up Mode, for lab use with automatic ramping up illumination to minimizes composite shrinkage on curing the indirect restorations e.g. crown, bridges, veneer, onlay or inlay.
- M3** : Pulse Mode, switching with high and low light output in less heat effect for clinical flexibility using for curing liner, bonding adhesives and the resin cement pre-cure.
- : Blue light and UV light
- : Full Battery
- : Battery Charging
- : Overheating - Device will automatically shut off and will be able to restart after cooling down.

Charger Signals:

Place the charger stand on a flat stable surface. Connect the power cable to the charger stand and plug into an electricity outlet, the green ready light will light up.

Note: If the battery capacity is considerably reduced and cannot be fully charged even after a long time, the battery is depleted. If the battery is depleted or damaged, please contact your dealer to buy a new battery.

CAUTION: Lithium Battery

1. Please charge and store the battery in a cool environment.
2. If the battery is warm after charging, wait for the battery to cool before using to maximize the battery life.
3. Avoid storing the battery in a fully charged state, as it will shorten the battery lifespan. If the battery is fully charged, emit the unit for about 20 seconds before storing.
4. Avoid storing the battery when nearly empty as self-consumption will result in excessive discharge and prevent further charging. If the battery is empty, charge for 10-20 minutes before storing.
5. When storing the battery, do not place on the charger stand to avoid overcharge.
6. Please charge/discharge the battery every 2 months when not used.
7. Please discard the battery according to manufacturer's instructions and recycle the battery.
8. Do not dispose in fire to prevent injury or burns.
9. Avoid the battery from contacting any metal objects or it will short circuit. Remove and reconnect the battery module if any short circuiting occurs.
10. Do not disassemble.
11. Do not crush.
12. Keep the battery dry.
13. Keep the battery away from the heat source.

Troubleshooting:

DESCRIPTIONS		
FAILURE	CAUSE	SOLUTION
Device doesn't power on	The battery is not charged	Check if the battery has been fully charged
	Defective battery or power PC board	Return to BA/ approved dealer for repair
LED lamp is unable to emit	Defective LED module	Return to BA/ approved dealer for repair
Handpiece doesn't charge	Poor connection between the charger stand and power adapter	Check the connection between the charger stand and power adapter
	The handpiece has been placed in the charger stand incorrectly	See the user manual for correct placement
	There may be moisture or debris in the charger stand	Clean out any liquid or debris from the charger stand
	Defective charger stand or adapter	Return to BA/ approved dealer for repair
Handpiece display panel doesn't light up	The battery is not charged	Check if the battery has been fully charged
	Defective LCD PC board, power PC board, or battery	Return to BA/ approved dealer for repair
Charger stand display panel doesn't light up	Poor connection between the charger stand and power adapter	Check the connection between the charger stand and power adapter
	Defective PC board, charging board, or adapter	Return to BA/ approved dealer for repair

※ When requested to send the unit back for repair, please send all the accessories (handpiece, charger stand, and power adapter) for further inspection.

Specifications:

Model number	BASE280
LED lamp	5 W
Light intensity	M1: 2,000 mW/cm ² and up M2: 1,600 mW/cm ² and up M3: 1,600 mW/cm ² (strong intensity) and up 600 mW/cm ² (weak intensity) and up
Wavelength range	Blue Light: 420 ~ 490 nm (peaks at 455 ~ 465 nm) UV Light: 360 ~ 420 nm (peaks at 400 ~ 410 nm)
Handpiece dimensions	250(L) x 40(W) x 90(H) mm
Handpiece weight	215 g
Fiber-optic probe dimensions	75 mm(L) x Ø10 mm
Fiber-optic probe weight	25 g
Battery specifications	7.4V, 1350 mAh
Charger stand dimensions	150(L) x 90(W) x 120(H) mm
Charger stand weight	260 g
Adapter power input	AC100~240V 50/60Hz
Adapter power output	DC 9V, 1.3A
System operational conditions	Temperature: 25°C to 35°C (77°F~95°F) Humidity: 15% to 95% @ 40°C, noncondensing Atmospheric Pressure: 700-1060 hpa
System storage conditions	Temperature: -20°C to 65°C (-4°F~149°F) Humidity: 15% to 95% @ 40°C, noncondensing Atmospheric Pressure: 700-1060 hpa
System transportation conditions	Temperature: -20°C to 65°C (-4°F~149°F) Humidity: 15% to 95% @ 40°C, noncondensing Atmospheric Pressure: 700-1060 hpa

EMC Tables:

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The model BASE280 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model BASE280 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Recommended separation distances between portable and mobile RF communications equipment and the model BASE280			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity			
The model BASE280 is intended for use in the electromagnetic environment specified below. The customer or the user of the model BASE280 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95% dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95% dip in UT) for 5 sec	<5 % UT (>95% dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model BASE280 requires continued operation during power mains interruptions, it is recommended that the model BASE280 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The model BASE280 is intended for use in the electromagnetic environment specified below. The customer or the user of the model BASE280 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BASE280, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$ <p>.....</p> <p>$d = 1,2 \sqrt{P}$ 80 MHZ to 800 MHZ</p> <p>.....</p> <p>$d = 2,3 \sqrt{P}$ 800 MHZ to 2,5 GHZ</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model BASE280 is used exceeds the applicable RF compliance level above, the model BASE280 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model BASE280.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than V/m.</p>			

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