
Maintenance and infection control

After every patient

Surfaces

- **1** Wipe with Planmeca approved surface disinfectant:
 - · metal parts of headrest
 - · metal parts of armrests
 - protective fascia of the operating light
 - handles of the operating light
 - hygienic membrane
 - instrument hoses
 - · instruments (see instrument instructions)
 - instrument console*
 - balanced instrument arms
 - tray table
 - · cuspidor top cover
 - · outside of the bowl
- Wipe the upholstery with mild soap- and water solution. Wipe with Planmeca approved upholstery disinfectant if needed.
- **3** Wipe the monitor and the handle (see monitor instructions).
- Oclean the bowl:
 - Pour a few drops of Planmeca approved suction disinfectant into the bowl and use a soft brush.
 - · Rinse the bowl by pressing the bowl rinse button.




*Before cleaning the instrument console, lock the touch screen.

- 6 Touch screen is unlocked. Press to lock.
- Touch screen is locked. Press and hold to unlock.

Before covering the instrument console and touch screen with a protective cover (optional), make sure they are dry.

Suction system


- Remove used aspirating tips.
- 2 Rinse the suction tubes by aspirating one glass of water with each suction tube.
- **3** Wipe the suction handpieces with Planmeca approved surface disinfectant.
- Wipe the suction holder with Planmeca approved surface disinfectant.



Waterlines

Short flushing

- Place all water consuming instruments into the flushing holder.
- 2 Bend the instrument arms to an angle of at least 90° and start short flushing.



In the morning

Suction system


- Rinse the suction tubes by aspirating 0.5 I of water with each suction tube.
- 2 Wipe the suction handpieces with Planmeca approved surface disinfectant.





Filters


Clean the coarse filters.



Waterlines

Long flushing

- Place all water consuming instruments into the flushing holder.
- 2 Bend the instrument arms to an angle of at least 90° and start long flushing.






After working day

Surfaces


Wipe the upholstery with Planmeca approved upholstery disinfectant.



Suction system

Suction cleaning

- Insert the suction tubes into the suction tube cleaning holder.
- 2 Start suction cleaning.
- 3 Disinfect the suction handpieces in a washer-disinfector.
- Suction holder: Disinfect the individual holders in a washer-disinfector. Wipe the intraoral scanner holder with Planmeca approved surface disinfectant.




Bowl

• Empty, clean and disinfect the bowl filter in a washer-disinfector.

Attention! Do not empty the bowl filter to the drain.


2 Remove the bowl and wash at max 65°C.



Waterlines

Long flushing

- Place all water consuming instruments into the flushing holder.
- Bend the instrument arms to an angle of at least 90° and start long flushing.




3 Switch off the unit.

Weekly

Surfaces

- Wipe all surfaces with mild soap- and water solution.
- 2 Wipe the upholstery with Planmeca approved upholstery treating agent.



Suction system

Dürr separators or VS/A systems only

- Mix 5 ml of Planmeca approved suction disinfectant with 250 ml of water.
- 2 Press the programming button for 5 s to open the bowl valve for 15 s to remove excess water.
- **3** When the valve has closed, pour the mixture into the bowl.
- 4 Let the solution affect overnight.





Waterlines

Waterline cleaning

- Place all water consuming instruments into the flushing holder.
- 2 Bend the instrument arms to an angle of at least 90° and start waterline cleaning.
- 3 Close the water tap.
- Fill the container with Planmeca approved waterline disinfectant. Attach the container to the unit. Follow the instructions on the display.
- **3** The next morning, switch the unit on.
- Follow the instructions on the display.

Chlorine check for Planmeca ActiveAqua™

- Measure the concentration of free chlorine in the water when H90.1 is displayed.
- 2 Select the corresponding value on the control panel.



Monthly

Suction system and waterlines

- Disinfect the instrument flushing holder in a washer-disinfector. Autoclave if needed.
- 2 Disinfect the suction tube cleaning holder in a washer-disinfector.
- 3 Empty the deposit cup when 70% full (WET systems only).

1 2 3

Filters

- Replace the coarse filters.
- 2 Empty the VS/A fine filter.
- 3 Empty the oil collector.

