
### CERASMART<sup>™</sup>

# The new leader in hybrid ceramic blocks



Introducing CERASMART<sup>™</sup>, a new force absorbing hybrid ceramic block




Contains high density of ultra fine homogeneously-dispersed fillers in a highly cross-linked resin matrix.


The genius behind GC's force absorbing hybrid ceramic is the development of a unique technology to manufacture a block with high density of ultrafine glass particles. Uniform dispersion (very short inter-particle distance) of individually silanated and bonded particles is key to delivering CERASMART's<sup>™</sup> exceptional strength, polish retention and wear resistance.



Dr Anthony Mak (Australia), Brad Grobler (New Zealand)



Typically contains nano ceramic particles embedded in a highly cross-linked resin matrix.




Typically contains a ceramic network strengthened by a polymer network.

# CERASMART<sup>™</sup> offers your patients strong, aesthetic and highly wear-resistant restorations

#### Perfect margins

CERASMART's<sup>™</sup> sharp margin adaptation and resistance to chipping ensures a long-lasting marginal seal.





CERASMART<sup>™</sup> Feldspathic Ceramic

Lithium Disilicate Glass-Ceramic


cate Hybrid

Hybrid Ceramic

**Resin Nano Ceramic** 

### High fracture toughness


CERASMART<sup>™</sup> has a very high level of breaking energy.



#Breaking energy is a measure of fracture toughness. R&D Dept. GC Corporation.

### Superior wear resistance

Independent research showed that  $CERASMART^{m}$  has very low wear loss under laboratory testing.<sup>\*</sup>



 \* Pande R et al. In Vitro Abrasion Wear of Nano and Hybrid Ceramic CAD/CAM Materials. J Dent Res 2015; 94(Spec Issue A):3529.
\*\*Not trademarks of GC Corporation.

# CERASMART<sup>™</sup> is easy to finish and retains its beautiful high gloss over time

Due to the dynamic proprietary nano ceramic matrix, and its complete homogeneous nature, CERASMART<sup>™</sup> is a true selfpolishing material. Not only does it stay polished longer, but it is proven to gain lustre even after being roughened.<sup>†</sup> This unique hybrid ceramic is unsurpassed with its high gloss value, unique selfpolishing capability, and its unmatched aesthetics.



Milled CERASMART™




CERASMART<sup>™</sup> after polishing



CERASMART<sup>™</sup> finished with OPTIGLAZE<sup>™</sup> color



Dr Anthony Mak (Australia), Brad Grobler (New Zealand)



+ Akiyama S, Akatsuka R, Sasaki K. Evaluation of Gloss-retention and Self-polishing Property of CAD/CAM Composite Block. J Dent Res 2015;94(Spec Issue A):3650.

### Finishing CERASMART<sup>™</sup>

### Try in and finish...



Remove the sprue of the restoration.



Try-in the restoration and check the occlusion.



Finish and polish using medium and fine silicone points.



Final polish using diamond polishing paste, eg. DIAPOLISHER PASTE.

### ...or characterise with OPTIGLAZE<sup>™</sup> color

### 1. Prepare



Sandblast with 25-50µm alumina (0.2MPa).

#### 2. Clean



Clean with oil-free air syringe or ultrasonic cleaner. Clean further with alcohol.

### 3. Silanate



Apply silane coupling agent, eg. CERAMIC PRIMER II and dry.

#### 4. Characterise



Dispense and brush apply OPTIGLAZE<sup>™</sup> color in a thin layer. Do not air blow. MUST be cured with a light-curing device with a wavelength in the range of 400-430nm (NOT with traditional Blue LED).

### Bonding CERASMART<sup>™</sup> Essential bonding steps for optimum results

### 1. Sandblast or 5% HFl etch



Options for pretreating CERASMART<sup>™</sup>: A. Sandblast with 25-50µm alumina (0.2MPa); OR

B. Treat with 5% HFI (hydrofluoric acid) for 60

seconds. Wash with water and dry.

2. Clean



Clean with oil-free air syringe or ultrasonic cleaner and dry. Clean further with ethanol (alcohol) to remove oil residue.

### 3. Silane




Apply silane coupling agent, eg. CERAMIC PRIMER II and dry. No bonding agent is required after the application of CERAMIC PRIMER II. For other ceramic primers, follow the manufacturer's instructions.

## Key elements for durable bonding

## to CERASMART<sup>™</sup>:

- 1. Surface roughness
- 2. A clean bonding surface
- 3. Applying CERAMIC PRIMER II

### Bonding to CERASMART<sup>™</sup> – effect of pre-treatment and thermocycling (TC)



When bonding G-CEM LinkAce resin cement to CERASMART<sup>™</sup>, similar bond strengths are achieved by either treating CERASMART<sup>™</sup> with 5% Hydrofluoric acid for 60 sec. and CERAMIC PRIMER II or by sandblasting and applying CERAMIC PRIMER II.

### **Completing Cementation**

**CERAMIC PRIMER II** features GC's advanced silane and proven phosphoric acid ester monomer technologies in an innovative single bottle delivery, formulated for stability and durable adhesion. It has a 2 year shelf life and does not require refrigeration.

**CERAMIC PRIMER II** is designed to create a strong bond between all aesthetic indirect restorations and resin-based cements.









### 4. Tooth preparation

Clean and treat the prepared tooth surface following the instructions of the manufacturer of the adhesive resin cement.

Retaining a wide band of enamel, which is etched for mechanical retention, is ideal for optimum adhesion.

### 5. Adhesive resin cementation

Follow the procedure described by the manufacturer of the adhesive resin cement.

**6. Finish** Finish and polish the margins

### CERASMART<sup>™</sup> clinical cases





Dr G Koike, Japan, CERASMART<sup>™</sup> A1 HT





Dr R Rosenblatt, USA, CERASMART™ A2 LT

### $\mathsf{Q}$ &A

### 1. Why is CERASMART<sup>™</sup> so aesthetic?

CERASMART<sup>™</sup> ultrafine glass filler particles allow the transmission of light in a way that reveals the natural opalescence, ensuring natural shade integration for highly aesthetic restorations.

#### **Opalescence (Transmitted light)**



Resin Nano Ceramic

Hybrid Ceramic

#### Fluorescence (Black light)






Resin Nano Ceramic Hybrid Ceramic


### 2. How easy is it to achieve and maintain a high gloss on CERASMART<sup>™</sup>?

CERASMART<sup>™</sup> is very easy to polish and due to its homogeneous, highly dense matrix, CERASMART™ will self-polish and gain lustre over time.

#### Gloss change with tooth brush wear test (1 year simulation)\*



#### 3. How radiopaque is CERASMART<sup>™</sup>?



\*R&D Dept. GC Corporation

### Q&A (continued)

### 4. What shades and sizes are available?


#### High Translucency (HT) Shades

Mainly used to replace enamel; use in cases of inlays, onlays, veneers, partial and full crowns



#### Low Translucency (LT) Shades


For replacement of dentin and enamel structures; use to mask discoloured preparations, in particular for crowns



#### Sizes, dimensions (mm) and indications Width Size Length Height Indications Size 12 15 12 10 More suited for inlays, onlays 18 12 More adapted for Size 14 14 individual crowns More suited for larger Size 14L 18 14 14 restorations (endocrowns, canines)

### 5. What is the flexural strength of CERASMART<sup>™</sup>?


CERASMART<sup>™</sup> features the highest flexural strength in its category, and is superior to the classical feldspathic ceramic blocks. Most importantly, it also offers a high breaking energy to buffer masticatory pressure and sustain more challenging clinical situations.



R&D Dept. GC Corporation

### 6. What are the preparation guidelines for CERASMART<sup>™</sup>?

- Margin preparation: deep chamfer or rounded shoulder.
- Prepare tooth with about 6° taper
- All internal edges and angles should be rounded
- Avoid having margins in direct occlusal contact with the opposing tooth



### 7. How can I stain and characterise CERASMART<sup>™</sup> restorations?

For external staining of CERASMART<sup>™</sup> restorations, use OPTIGLAZE<sup>™</sup> color. OPTIGLAZE<sup>™</sup> color is a light-cured, nano-filled characterisation system specially designed to add colour, surface gloss and high wear resistance to CERASMART<sup>™</sup> restorations.

Once applied, OPTIGLAZE<sup>™</sup> color's single dispersion nano-filler technology ensures a toughened 25-50 µm characterisation layer to give your restorations beautiful, colour-stable aesthetics and a long lasting gloss with exceptional wear resistance.

#### Nano-filler technology





### 8. What polishing paste can I use with CERASMART<sup>™</sup>?

DIAPOLISHER PASTE is GC's polishing paste, featuring ultrafine diamond particles (1µm). It is used as a final step after your usual polishing procedures and will provide a very high gloss on any direct or indirect composite restoration.



### 9. How do I repair CERASMART<sup>™</sup> restorations?

- 1. Slightly roughen the bonding area (restoration and exposed tooth structure) using a coarse diamond point or carbide bur.
- 2. Apply silane coupling agent on the restoration such as CERAMIC PRIMER II and gently dry. If another primer is used, follow the manufacturer's instructions.
- 3. Apply a bonding agent, eg. G-ænial Bond, to all tooth and restoration surfaces to be repaired. Dry and light-cure following manufacturer's instructions.
- 4. Apply the chosen restorative material and light-cure according to the manufacturer's instructions.

### Product Range



#### CERASMART<sup>™</sup>

CERASMART<sup>™</sup> for CEREC, refill of 5 blocks Sizes 12, 14, 14L HT shades (A1, A2, A3, A3.5, B1) LT shades (A1, A2, A3, A3.5, B1) Bleach shade (BL)



CERAMIC PRIMER II 3ml bottle



OPTIGLAZE<sup>™</sup> color Complete set



DIAPOLISHER PASTE 2g syringe



GC Australasia Dental Pty Ltd 1753 Botany Road, Banksmeadow NSW 2019, Australia T: +61 2 9316 4499 F: +61 2 9316 4196 www.gcasia.info