


UNLIMITED

coltène®

Lab-Putty







Key for gingival mask



Model for repairs

	1	
Indications	Lab-Putty	Lab-Putty hard
Matrix for the formation of general occlusion and functional patterns	()	Y
Model for repairs	()	()
• Soldering		(
Isolation / denture heat protections	()	
Hard palate relief stamps	0	0
• Key for the gingival mask	0	
Block outs	•	
Matrix to built partial and full dentures	7 O 🔀	
• Relinings	71	0

TECHNICAL DATAS LAB-PUTTY: *

Product	Mixing time in s	Working time in s	Setting time in s		Dimen stabili	A SECTION ASSESSMENT	Heat resistence in °C
Lab-Putty	30-45	90	5-7	63-70	1h 1day	- 0.08 - 0.16	110 °C 230 °F

^{*} after dosage according to instructions for use

easy and safe dispensing / non-irritating paste activator / convenient mixing due to its consistency / flexible working due to t





coltène®

Lab-Putty hard







Matrix for facings



Matrix for full dentures

Heat resistant up to 300 °C



Broken partial denture



Soldering of the point of fracture



After polishing

Coltène Laboratory modeling materials: Lab-Putty and Lab-Putty hard used in the appropriate applications have important benefits for the lab-technician:

Easy and safe dosage, convenient mixing consistency, problem-free removal and an ideal elasticity after setting. In addition, the optimised end hardness and out-standing dimensional stability provide perfect accuracy and excellent results in a time saving manner.

During the use of both types of Lab-Putty for denture relining, keys and repair matrices, the step of separate isolation is eliminated and a perfect smooth acrylic denture surface is guaranteed.

Testimonial: «ZahnArt» dental-laboratory Dornbirn, Austria

TECHNICAL DATAS LAB-PUTTY HARD: *

Product	Mixing time in s	Working time in s	Setting time in s	Endhard- ness in shore A	Dimension stability in %	Heat resistence in °C
Lab-Putty hard	45	120	15	90	1h - 0.04 1day - 0.13	300 °C 570 °F

* after dosage according to instructions for use

he unique activator / problem-free removal as a result of ideal end hardness and excellent elasticity / easy and safe dispensing /





coltène®

Lab-Putty

coltène®

Lab-Putty hard



Order information

Lab-Putty

8805 Lab-Putty Base Activator 40 ml Dosing spoon 1

8815 Lab-Putty Base 7.5 kg / 5000 ml Dosing spoons 5

8840 Lab-Putty / Lab-Putty hard Activator

Lab-Putty hard

8882 Lab-Putty hard Base 1.800 kg / 900 ml Activator 40 ml Dosing spoon 1

8880 Lab-Putty hard Base 5 kg / 2600 ml Dosing spoons 5

8840 Lab-Putty / Lab-Putty hard

Activator 40 ml



Coltène/Whaledent Ltd.

40 ml

The President Suite Kendal House, Victoria Way Burgess Hill, West Sussex RH15 9NF / UK Tel +44 1444 235486 Fax +44 1444 870640 info.uk@coltene.com

Coltène/Whaledent AG

Feldwiesenstrasse 20 9450 Altstätten/Switzerland Tel. +41 71 757 53 00 Fax +41 71 757 53 01 info.ch@coltene.com

Coltène/Whaledent Inc.

235 Ascot Parkway Cuyahoga Falls, Ohio 44223/USA Tel. +1 330 916 8800 Fax +1 330 916 7077 info.us@coltene.com

