TISSUE MANAGEMENT



* Data on file. ** Trademark of a company other than Ultradent. 1. realityesthetics.com. 2. "Can't Live Without" Clinical Research Associates Newsletter, Volume 21, Issue 7, July 1997.

aluminum chloride, etc. are also enhanced when using Ultrapak plain knitted cords, which carry significantly greater quantities of hemostatic

solution than conventional cords.

TISSUE MANAGEMENT

PRE-PREPARATION PACKING TECHNIQUE

To ensure cord retention during preparation, use a cord large enough to firmly compress into sulcus.

PREPACK



1. Place Ultrapak[™] knitted cord soaked in hemostatic solution using a cord size that appears slightly too large to ensure cord retention. The thin Ultrapak[™] Packer quickly slips cord into position. The knitted cord's unique design (interlocking loops) facilitates easy packing and locks it into place.

PREPARATION



2. Extend margin subgingivally by cutting partway into knitted cord, which won't entangle in diamond bur. Remove remnant of cord with an explorer or other instrument. Bleeding is minimal if at all. A small portion of uncut tooth above gingival attachment is preserved to record in impression. If additional retraction is required, repack with appropriately sized cord. Rinse, air dry, and make impression.

DOUBLE-CORD TECHNIOUE

The most common challenges in getting a quality impression are adequate tissue retraction and sufficient moisture control. Try a double retraction cord technique combined with effective hemostatic agents to alleviate both.

FIRST CORD



1. Once hemostasis is achieved, carefully place a single cord —such as Ultrapak" knitted cord #0, #00, or #000—to the bottom of the sulcus. Use Fischer's Ultrapak Packers to place cords properly and efficiently.

SECOND CORD



2. Place a second, thicker cord soaked in a hemostatic agent to expand the tissue laterally.

RINSE/DRY



3. Rinse the area well, lightly dry, and make impression.





Knitted Ultrapak cord is composed of thousands of tiny, interlocking loops so it compresses and expands easier than other cords. 100% cotton fibers provide high absorption of hemostatic agents and sulcular fluids.



Ultrapak cord, saturated with hemostatic solution, controls bleeding and sulcular fluid near gingival and subgingival preparations.

0800 808 855



Ultrapak cord compresses upon packing and then expands for optimal tissue displacement.



Ultrapak CleanCut design features a blade in the cap for efficient cutting. A special dispensing orifice prevents cord from falling into the bottle.

FOR DIGITAL IMPRESSIONS -COMPLETE HEMOSTASIS

HEMOSTASIS



1. Complete hemostasis is essential, especially when taking digital impressions, for the most accurate marginal fit of any restoration.

CLEAR FIELD



2. After hemostasis is achieved and tissue is retracted, preparation is ready for digital impression.

"Ultrapak cord is excellent at displacing the gingival tissue and allowing proper hemostasis, and is easy to place and remove. It works very well for all retraction purposes. The different sizes are good for all situations." —DR. Y ČLEMĚNT SHEK – SAN FRANCISCO, CA

"In dentistry, time is money. Ultrapak cord's woven design makes packing the cord quick and easy, plus the tooth can be prepped or touched up without snagging the cord. This increases patient comfort in shortening the appointment with far less repeat impressions." —DR. THOMAS J. FRANKFURTH – TAMPA, FL

"Ultrapak cord has taken the stress out of cord packing. This was the most frustrating part of my day when I was using other products. Add the amazing Astringedent hemostatic and... LIFE IS GOOD!"

—DR. LISA MARSHALL – XENIA, OH

