



## A New Soft Tissue Retractor— So Comfy Your Lips Will Thank You

### STEP 1.



Simply press the tabs on the Umbrella™ soft tissue retractor together to prepare it for insertion.

### STEP 2.



Ask the patient to place the tip of their tongue on the roof of their mouth.

### STEP 3.



Choose one side of the mouth in which to start, and then comfortably insert the other side of the retractor into the cheek.

### STEP 4.



Use the tabs to center the retractor with the patient's mouth.

### STEP 5.



Check that the patient's tongue is resting comfortably behind the guard, ensuring easy access.

The new Umbrella retractor has been designed for your patient's comfort, so that you can focus all attention on the procedure at hand. It couldn't be easier to place. Just slide under the lips, and the tongue sits gently behind the built-in rest. Try it yourself, and just for fun, see how well you can say these phrases:



- ◆ Umbrellas Redirect Raindrops
- ◆ Rest, Retraction, and Rubber Baby Buggy Bumpers
- ◆ Soft on Your Cheeky and Sneaks in Sneaky
- ◆ Tongue Talk Used to Be Tough
- ◆ Lips Love to Let Loose