PHILIPS

Sonic electric toothbrush with app

DiamondClean 9000

Built-in pressure sensor Smart brush head recognition 4 modes, 3 intensities Connected brushing made easy





Whiter, healthier teeth for life Whiter teeth in just 1 day*

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

Enjoy a stylish and user-friendly design

• Simple to charge and style to boot

DiamondClean users experience improved oral health

- Hard on plaque, gentle on your gums
- Removes up to 10x more plaque* for a deep clean

Let the DiamondClean 9000 adapt to your needs

- Personalized brushing experience
- · Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- BrushSync automatically selects the best mode for you
- · Always know when to replace your brush heads

Stay on track with a personalized progress report

• Start and keep up healthy habits

Highlights

Removes up to 10x more plaque*



Know you're getting the deepest possible clean with The C3 Premium Plaque Control brush head. The soft, flexible bristles are designed to curve around the contours of your teeth, giving you 4x more surface contact^{**} and up to 10x more plaque removal from hard-toreach spots.

C3 Premium Plaque Control



The C3 Premium Plaque Control brush head is designed to give you your deepest clean yet. Soft flexible sides and bristles perfectly curve around the contours of your teeth, giving you 4x more surface contact and helping target hard-to-reach spots^{*}.

Progress Report



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

4 modes, 3 intensities



DiamondClean 9000 comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

Long-life battery



Your DiamondClean 9000 will last for up to two weeks on a single charge with regular use.

BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

Specifications

Modes

3 intensities: High, Low, Medium Clean: For exceptional everyday clean Deep Clean+: For an invigorating deep clean Gum Health: Special attention to molars White+: To remove surface stains

Items included

Handle: 1 DiamondClean 9000 Brush heads: 1 C3 Premium Plaque Control, 1 W2 Optimal White Glass charger: 1 Travel case: USB charging travel case Charger base: 1

Connectivity

Bluetooth® wireless technology: Connected brushing app

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

Design and finishing

Color: White

Cleaning performance

Performance: Removes up to 10x more plaque*

Whitening benefits: Whitens teeth in 1 day* Pressure feedback: Vibrates handle to alert user

Timer: BrushPacer and SmarTimer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard BrushSync Replacement Reminder: Always know when to, replace brush heads BrushSync Technology: smart brush head, Connects smart handle and

Ease of use

Handle compatibility: Easy click-on brush heads Replacement reminder: To always ensure best results, reminder icon lights up Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable Battery type: Lithium ION Operating time (full to empty): 14 days***

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

© 2020 Koninklijke Philips N.V. All Rights reserved.

their respective owners.

Issue date 2020-08-03 Version: 4.0.1

EAN: 08 71010 39326 59



* * compared to DiamondClean

* ** based on two periods of two-minute brushings per day

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or

www.philips.com

^{*} more than a manual toothbrush