




Philips Sonicare  
ProtectiveClean 6100  
Sonic electric toothbrush

Built-in pressure sensor

3 modes, 3 intensities  
2 x BrushSync features  
Travel case



HX6877/23

## Whiter teeth. Gently does it.

Whitens teeth in just 1 week.

Feel the difference of a gentle clean with our Pressure Sensor while whitening your teeth in 1 week.

### Proven to improve oral health

- Whitens teeth in just one week
- Safe & gentle on sensitive areas, orthodontics & dental work

### Optimize your brushing

- Three modes, three intensity settings
- Brush heads select the optimal mode and intensity

### Innovative technology

- Lets you know when you're pressing too hard
- Connects smart brush handle and smart brush heads
- Always know when to replace your brush heads
- Dynamic cleaning action drives fluid between teeth

### Designed around you

- Encouragement to brush thoroughly
- Makes traveling easier for you

**PHILIPS**  
**sonicare**

# Highlights

## Naturally whiter teeth



Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stain-removal bristles, it's clinically proven to whiten teeth in just one week.

## Brushing your way



This toothbrush allows you to customize your brushing with a choice of three modes and three intensities. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums. Three intensities allow you to switch the order of 'higher' and 'lower' to align with the intensities from left to right.

## Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

## Philips Sonicare technology



Your Sonicare brush head and handle are the ultimate team. Enhancing our sonic technology, they give you unrivaled cleaning. High-speed vibrations in your toothbrush handle power up to 62,000 brush movements per minute that whip up your toothpaste into gentle cleaning bubbles, driving them deep between your teeth and along the gum line for exceptional cleaning.

## Too much pressure?



You may not notice if you're brushing too hard, but your toothbrush will. If you start taking it out on your gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

## BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart mode pairing and Smart replacement reminders.

## Mode and intensity pairing



Wondering which mode and intensity you should use? Wonder no longer. BrushSync mode pairing tells your smart handle which smart brush head you're using. So, if you click on a gum care brush head, for example, your toothbrush knows to select the optimal mode and intensity for your gums. All you need to do is to press the power button.

## BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

## Helpful timers



Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

# Specifications

## Modes

- Clean: For exceptional everyday clean
- Gum Care: Gently massages gums
- White: Removes surface stains
- 3 intensities: Low, Medium, High

## Items included

- Handles: 1 ProtectiveClean
- Brush heads: 2 W2 Optimal White standard
- Travel case: 1
- Charger: 1

## BrushSync Mode Pairing

- C2 Optimal Plaque Control: Pairs with Clean mode
- G2 Optimal Gum Care: Pairs with Gum Care mode
- W2 Optimal White: Pairs with White mode

## Design and finishing

- Color: White Silver

## Cleaning performance

- Speed: Up to 62000 brush movement/min
- Health benefits: Improves gum health in only two weeks, Helps reduce cavities

- Pressure feedback: Vibrates handle to alert user
- Timer: Smartimer and Quadpacer

## Smart sensor technology

- Pressure sensor: Alerts when brushing too hard
- BrushSync Replacement Reminder: Always know when to, replace brush heads
- BrushSync Technology: Connects smart handle and, smart brush head

## Ease of use

- Handle compatibility: Easy click-on brush heads
- Battery indicator: Light shows battery status
- Handle: Slim ergonomic design
- Brushing time: Up to 2 weeks

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Operating time (full to empty): Up to 2 weeks

## Power

- Voltage: 110-220 V

## Service

- Warranty: 2-year limited warranty



Issue date 2018-06-27

Version: 3.1.1

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* than a manual toothbrush