

B Classic

Designed with performance, Reliability & Safety in mind.



- LCD display with keypad and test reminder. Never forget to run important Vacuum & Helix tests.
- Performance vacuum pump for speed and reliability
- Optional Barcode and cycle report printer. Track critical instruments digitally.
- Delayed start so Vacuum and Helix tests can be programmed to run before a practice opens. A real time saver.
- Wi-Fi and Ethernet port are standard inclusions so Steriliser can be networked.
- User management system allows users to be tracked to cycle data via personal PIN codes.
 Complies with Australian and NZ standards and means no paper printout is required.
- Cycle data saved as a PDF file meaning no software required to save data. Can be downloaded to PC via included USB.
- Steriliser can be installed as a stand alone device or connected to external R/O systems & Autofill options. Multiple solutions to cover varying installation requirements.
- Built in water conductivity sensor to ensure water quality is always at it's best.
- Software available to manage automatic upload of cycle data. Traceability software available to track critical instruments to the patient.

Industry leading 3 year or 3,000 cycle warranty!