

EDGE ONE

FIRE™

THE NEW STANDARD...
ONE THIRD THE COST.



\$63.00
per 3 pack

EdgeOne Fire™ is designed to shape canals in a reciprocating motion similar to WaveOne® Gold and is available in sizes Small, Primary, Medium and Large. EdgeOne Fire™ can be used with the same motor and hand piece setting as WaveOne® Gold. EdgeOne Fire™ features our heat-treated Fire-Wire™ NiTi and tests at five times the cyclic fatigue as WaveOne® Gold. Experience the new standard - at one third the cost!

Tip Size	21mm	25mm	31mm
20	EN-EOFSM21	EN-EOFSM25	EN-EOFSM31
25	EN-EOFPRI21	EN-EOFPRI25	EN-EOFPRI31
35	EN-EOFMED21	EN-EOFMED25	EN-EOFMED31
45	EN-EOFLG21	EN-EOFLG25	EN-EOFLG31



EDGEENDO

HENRY SCHEIN®
DENTAL

Rely on Us™