

POWERFUL YET GENTLE: LUNOS® PASTES



LUNOS® PROPHY PASTE TWO IN ONE

The innovative two-in-one solution: removes discolouration and polishes in a single step. Saves time, yet is thorough and gentle on tooth substance.

The abrasiveness of the paste decreases gradually thanks to the self-reducing polishing particles

- Can be used for all treatment needs
- Effective cleaning and polishing thanks to an innovative system of abrasive and polishing agents. Hydroxyapatite supports remineralisation
- Particularly gentle on tooth substance, RDA 28
- Mint flavour additionally contains sodium fluoride to support remineralisation

Flavours: Mint. 100g tube



LUNOS® PROPHY PASTE SUPER SOFT

Lasting, gentle fine polish, even on particularly sensitive surfaces.

- Universal use, even on implants, orthodontic appliances and dental restorations
- Creates smooth, high-gloss surfaces, and new layers of plaque or calculus take longer to form
- Hydroxyapatite cleans and polishes extremely gently and supports remineralisation
- RDA < 5 "Relative Dentin Abrasivity"

Flavours: Neutral 50g tube



LUNOS® PROPHYLAXIS RING

For optimised handling through simple application of the propy paste.

- Flexible ring size
- Re-usable