

Surgic Prot × Vario Surg 3

Dual Power and Performance

NSK provides **the world's first solution** linking two essential elements of implant procedures—an implant motor and ultrasonic bone surgery device.

Surgic Pro, the implant drilling device boasting superior and precise torque control via AHC (Advanced Handpiece Calibration system)—and VarioSurg3, the new ultrasonic bone surgery unit offering 50 percent greater power—can link up and operate as a single device.

The dedicated Link Stand efficiently mounts two compact and stylish control units and requires only minimal space. The synergy effect of this will significantly expand oral surgery capabilities.







Surgic Pro+

- Powerful torque (up to 80Ncm)
- Wide speed range (200 to 40,000min-1)
- High torque accuracy
- Light and small micromotor with LED illumination
- Excellent durability
- Autoclavable and thermo-dininfectable micro motor

VarioSurg3

- Major power boost Up 150 percent
- The selectable "Burst" function generates a hammer effect
- \bullet Adjustable light intensity, water volume, and power mode.
- Slim and ergonomic handpiece with LED illumination
- Power feedback & auto tuning function
- A wide choice for clinical applications
- Autoclavable and thermo-dininfectable micro motor