

The leader of the PACK

MAXimize Your Vertical Compaction
with the **NEW Obtura MAXPACK**



Obtura Spartan Endodontics has offered dependable products to clinicians for over 25 years. The Obtura Maxpack is no exception.

The MaxPack offers:

- Lightweight, ergonomically designed handpiece allowing unlimited mobility
- Four temperature pre-sets ranging from 150°C to 230°C
- 12 month Warranty
- Obtura Spartan pluggers and spreaders designed to be used specifically with the Obtura MaxPack

To order your MaxPack TODAY, please call (800) 344-1321, visit us online at www.Obtura.com or email enhancingendo@obtura.com

 **OBTURA
SPARTAN**
Endodontics

MAXPACK Item #'s and Descriptions

824-700	Obtura MAXPACK *Includes MaxPack Unit, Charging Console, MaxPack Battery, Battery Enclosure, AC Power Adapter, Power Cord, and 1 Free Obtura Spartan 55/.06 Taper Plugger	824-702	MaxPack Battery
824-706	MaxPack Heating Unit	824-703	MaxPack Battery Enclosure
824-701	MaxPack Charging Console	824-704	MaxPack AC 12V 1620mAH Adapter
		824-705	MaxPack Power Cord

Obtura Spartan Pluggers

- Autoclaveable pluggers, made from tempered stainless steel
- Deliver consistent heat to soften gutta percha
- Available in ISO sizes and ISO tapers
- Also Available: Hot Testing Tip, Narrow/Standard Spreaders, and Tip Stand



Obtura Spartan Pluggers

● 824-710	Obtura Spartan 30/.04 taper Plugger
● 824-711	Obtura Spartan 35/.04 taper Plugger
● 824-712	Obtura Spartan 40/.04 taper Plugger
● 824-713	Obtura Spartan 45/.04 taper Plugger
● 824-714	Obtura Spartan 55/.06 taper Plugger
● 824-715	Obtura Spartan 55/.08 taper Plugger
● 824-716	Obtura Spartan 55/.10 taper Plugger
● 824-717	Obtura Spartan 60/.12 taper Plugger

Obtura Spartan Spreaders/Accessories

● 824-718	Obtura Spartan Hot Testing Tip
● 824-719	Obtura Spartan Spreader (Narrow) - 25/.03
● 824-720	Obtura Spartan Spreader (Standard) - 25/.05
824-722	Obtura Spartan Plugger Stand - Black (Holds 6)
824-723	Obtura Spartan Plugger Protector Sleeve (Bag of 10)

(800) 344-1321 | www.Obtura.com | enhancingendo@obtura.com

© 2010 Obtura Spartan