



STS Health

AquaClave MARK 3 RO P5

Why choose an AquaClave?

- **NO MORE STORING \$\$\$ BOTTLED WATER**
- **ENVIRONMENTALLY friendly**
- **CONSISTENT high quality water**
- **Automatic system ensures 10L is ready ON TAP**
- **ECONOMICAL choice for Perth water**
- **EASY to maintain, filters changed in-house**
- **AUSTRALIAN Designed Electronics**
- **USA Manufactured Housings**
- **Leak Detector - Easy to use & EFFECTIVE**

The AquaClave Mark 3 is an automatic multi-stage system using five purification steps to deliver exceptionally high quality water to your Sterilisation equipment.

The availability of high quality water on demand avoids the need for purchasing bottled water, avoids storage challenges, minimises steriliser maintenance costs and reduces the damage to instruments from tap water rinsing which makes buying an AquaClave an intelligent decision.



Reverse Osmosis System for Perth Water -

The AquaClave has been engineered to economically produce exceptionally high quality water even in the most challenging water quality environments around Australia.

Once plumbed the AquaClave automatically takes tap water and passes it through the first stage of advanced pre-treatment steps to remove any suspended solids, organic components and chlorine. Once complete the water is split through the Reverse Osmosis membrane. The remaining permeate continues on through a two stage polishing system to remove any remaining organic and non-organic residues thereby finishing the water to meet all autoclave operating standards.

Please note that this system is only suitable for use on disinfected potable water supplies from your local water authority. For areas where tank or other water is required please consult Aqualife. To achieve the water quality indicated above the system must be operated in accordance with the operating and maintenance instructions. On site verification of the water quality is required for municipal water supplies.

Supplied and Serviced by:



T: (08) 9244 4628

E: info@stshealth.com.au

www.stshealth.com.au

Manufactured by:



EXCELLENCE IN WATER TREATMENT
& LOW COST MAINTENANCE SYSTEMS

www.aqualife.com.au