


SciCan Sys™ Instrument Cassettes

Two Sizes Available

Sys™ is designed to fit in washers and autoclaves in a modular fashion, and ensures efficient use of space. Double hinged for easy lid management at chair side. Scalloped instrument supports keep instruments in place.



Sys-TM 1

- Ideal for exam set-ups
- Can hold up to 7 instruments
- 64 x 178 x 13 mm (2.5 x 7 x 0.5")



Sys-TM 3

- Ideal for perio/hygiene
- Capacity of 10 to 12 instruments
- 114 x 178 x 13 mm (4.5 x 7 x 0.5")

For efficient instrument cleaning, sterilisation and storage



Durable

High grade surgical stainless steel with a passivated finish to maximize durability and minimize discolouration.



Efficient

Efficient, time-saving transfer of instruments from the operator to HYDRIM® to your steriliser.



Warranty

Manufacturer's four-year warranty on all moving parts and a lifetime guarantee against corrosion.



Convenient

Keeps instruments organized and neat.