

# DRY MILLING WITH A PLUS.

## K5+



**5**

Axes

**1**

Disc

**6**

Blocks

**16**

Tools

**Dry**

**Mill**

**vhf**  
CREATING PERFECTION

# THE GLOBAL BESTSELLER. WITH EVEN MORE FEATURES.



## Unparalleled precision

- Restorations in ultra HD
- Premium spindle with 4-fold hybrid ceramic ball bearings for highest running accuracy
- 3 microns repetition accuracy



## Powerful and robust

- Mills the toughest materials on the market, incl. CoCr
- Powerful 500 watt spindle and 60,000 RPM
- Heavy industrial quality for maximum stiffness
- Solid cast body for the lowest vibrations



## Maximum independence

- Unlimited material accessibility in 98 mm disc format, separate block and prefab abutment holders available
- Covers the broadest range of indications, due to  $\pm 35^\circ$  rotation angle in the 5<sup>th</sup> axis, and up to 40 mm blanks



## Proven German reliability

- 100% engineered and manufactured in Germany
- Sophisticated sealing air concept to protect mechanics, electronics, and spindle
- Webcam for remote support
- Ethernet for stable long-range connection
- 24-month warranty



One-handed clamping made possible by **DirectDiscTechnology**



## Extremely economical

- Ionizer and improved air circulation for easy machine cleaning
- **DirectDiscTechnology** for revolutionary disc clamping
- Automatic changer for 16 tools
- Easy to operate DentalCAM with DirectMill function included – no annual fees

# TOOLS AND SOFTWARE


## vhf precision tools

- Grinding and milling tools with sophisticated cutting geometries
- For machining all material classes: wax, PMMA, zirconia, composites, titanium, cobalt-chrome, and glass-ceramics
- First-class service life, best quality and attractive prices

## vhf DentalCAM

- High-performance software package with open STL import for all common CAD programs
- Optimized CAM strategies for the fastest possible machining
- Perfectly adapted for all vhf dental milling machines – for predictable first-rate results

# PARTNERS WORLDWIDE



Dental milling machines, tools, and software from vhf are sold exclusively by official vhf resellers and OEM partners. Get in touch!

Global: [dental.vhf.de/en](https://dental.vhf.de/en)  
The Americas: [vhf.com](https://vhf.com)

■ Dental  
■ from  
■ Deutschland