## Introducing our <a href="new">new</a> Coxo Pens & Activators


## **Ultra Smart Al Endo Ultrasonic Activator**





38kHz frequency - Cleans the root canal more effectively

- In 3 Working Power mode choice



- Has good ultrasonic transmission function
- Autoclavable & in 3 sizes: 25/02, 30/04 & 30/00



Built-in ARM processor, can auto adjust the frequency according to the load. Detachable silicone cover

Safer working time limit - Auto stop in 20secs. 5sec beeps
Large Battery capacity - Sustainable standby 5 hours

Simple - One-key control of power on/off. Light = 650gm

## **C-Fill mini Endo Obturation Pen**



4 Sec Heating time

## **4 Precise working Temps**

150, 180, 200 & 230  $^{\circ}$ C For ALL Gutta Percha

Double Side Controls & Double Displays

