

BA Ulticlean Prophylaxis



# Powder series

**Ulticlean Prophylaxis  
Powder STANDARD**  
(pack of 4).  
For supragingival use.

**Ulticlean Prophylaxis  
Powder Perio PLUS**  
(pack of 4).  
For sub-and supragingival use.



# Powder information

|                                                                                     |                           | STANDARD           | STANDARD           | PERIO PLUS     |
|-------------------------------------------------------------------------------------|---------------------------|--------------------|--------------------|----------------|
| ITEM CODE                                                                           |                           | BA-LEMONPOWDER     | BA-MINTPOWDER      | BA-PERIOPOWDER |
| FLAVOUR                                                                             |                           | LEMON              | MINT               | NEUTRAL        |
| RAW MATERIAL BASE                                                                   |                           | SODIUM BICARBONATE | SODIUM BICARBONATE | TAGATOSE       |
| PARTICLE SIZE Ø                                                                     |                           | 40 µM              | 40 µM              | 15 µM          |
| pH                                                                                  |                           | AT 20 °C: 8,1-8,4  | AT 20 °C: 8,1-8,4  | AT 20 °C: 6-8  |
| TOOTH                                                                               | STAINS AND YOUNG TOOTH    | ✓                  | ✓                  | ✓              |
|                                                                                     | ENAMEL                    | ✓                  | ✓                  | ✓              |
|                                                                                     | WHITE SPOTS ON THE ENAMEL | ✓                  | ✓                  | ✓              |
|                                                                                     | DENTINE                   | ✓                  | ✓                  | ✓              |
| GUM                                                                                 | GUM                       | ✓                  | ✓                  | ✓              |
|                                                                                     | TONGUE AND PALATE         | ✓                  | ✓                  | ✓              |
| RESTORATION                                                                         | HYBRID COMPOSITE          | ✓                  | ✓                  | ✓              |
|                                                                                     | GLASS IONOMER             | ✓                  | ✓                  | ✓              |
| IMPLANT                                                                             | IMPLANT AND MUCOSITIS     | ✓                  | ✓                  | ✓              |
| ORTHODONTICS                                                                        | BRACKETS AND APPLIANCES   | ✓                  | ✓                  | ✓              |
|  |                           | ✓                  | ✓                  | ✓              |
|  |                           | ✓                  | ✓                  | ✓              |
|  |                           | ✓                  | ✓                  | ✓              |



243/249 Bush Road, Rosedale, Auckland 0632, New Zealand  
[www.henryschein.co.nz](http://www.henryschein.co.nz) | 0800 808 855



Scan &  
Order Online