



Shouldn't a tool you use every day be more comfortable?

No matter which member of the Filtek[™] family of products you're using, they can all be dispensed with the same device. That's why we've designed the 3M[™] Filtek[™] Restoratives Dispenser with improvements that can make your restorative delivery easier, more comfortable and secure.

Designed for warming

3M[™] Filtek[™] Restoratives Dispenser was designed with no metal parts, and can be used with compatible warming devices, with heating up to 70°C/158°F, for up to 1 hour. 3M has validated the following composites can be safely warmed: Capsules 3M[™] Filtek[™] Supreme XTE Universal Restorative, 3M[™] Filtek[™] Universal Restorative; Syringes 3M[™] Filtek[™] Supreme Flowable Restorative, 3M[™] Filtek[™] One Bulk Fill Restorative.



3M | Filtek

Easy to load

Secure capsule

Stop worrying about the

capsule shifting or falling out during use, thanks to a snap-in feature.

You can load the rotatable capsule holder at any angle, even using one hand.

Simple cleaning

You can clean it manually or with a washer-disinfector. It was safely tested for 200 cycles using high-level sterilisation and disinfection.

Extra comfort

Less straining on your hand, thanks to its lighter weight, easy-pull trigger and ergonomic design. It also reduces the force needed for extrusion.



3M[™] Filtek[™] Restoratives Dispenser works with ANY of 3M's composites.



3M™ Filtek™ Universal Restorative



3M™ Filtek™ Supreme XTE Universal Restorative



3M™ Filtek™ One Bulk Fill Restorative



3M™ Filtek™ Flowables



3M[™] Paradigm[™] Nano Hybrid Universal Restorative



3M™ Filtek™ Z250 Universal Restorative

3M.com.au/Filtek



3M Oral Care 3M Australia Pty Limited Building A, 1 Rivett Road, North Ryde NSW 2113 1300 363 878 3M.com.au/Dental

3M New Zealand Limited 94 Apollo Drive, Rosedale Auckland 0632 0800 808 182 3M.co.nz/Dental