


the all NEW

STRATA-G™

Now faster and more predictable!



Reliable, tight contacts

Drawn-wire nickel titanium plus PEEK produces consistent separating pressure longer!



Easier, more predictable contours

Marginal ridge enhancements built into the matrix bands AND rings help provide that finishing flair!

Faster, less frustrating Class IIs

The Strata-G™ rings, bands and wedges combine for a 270° Seal™ for ultimate flash control. The rings won't spring off the tooth, and the amazing Wide Preparation ring makes tough wide preps no problem.


270° Seal™ for ultimate flash control!



Garrison® sectional matrix systems are the standard of care.

Tofflemire-style System

- Fails to restore proximal anatomy
- Thin contact at the marginal ridge
- Large food trap (below)
- Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Strata-G™ Sectional Matrix System

- + Operator-friendly retaining system
- + Naturally contoured bands
- + Anatomically correct contacts
- + Contacts at the height of contour
- + Rings produce optimal tooth separation for excellent, tight contacts



faster—more predictable—flash control

STRATA-G™

Garrison® sectional matrix systems are the standard of care.

Longer lasting by design.

Advanced materials and design combine to produce our longest lasting, most predictable separator ring ever!



GA1-SG400
2 pk.



GA1-SG500
2 pk.



GA1-SG600
2 pk.

Anatomy top to bottom

Five sizes of super-thin stainless steel in both dead-soft and firm ensure you have the best matrix band to restore every case.

100 ct. bulk pack



SG050
3.5
100 pk

SG100
4.5 mm
100 pk

SG175
5.5 mm
100 pk

SG200
6.5 mm
100 pk


SG300
7.5 mm
100 pk

Also available in 50ct. bulk pack

Flexibility for a better cervical seal!

Strata-G™ Wedges utilize Soft-Face™ materials and advanced design to promote a better fit.

100 ct. bulk pack
Also available in 50ct. bulk pack



GA1-SG-KSH-10

Contents:

- 1 - Short Ring (blue)
- 1 - Tall Ring (orange)
- 1 - Wide Prep Ring (green)
- 50 - Assorted Matrix Bands
- 50 - Assorted Firm Matrix Bands
- 10 - VariStrips
- 80 - Assorted Ultra-Adaptive Wedges
- 5 - Deep Margin Elevation Bands
- 1 - Ring Placement Forceps

Garrison

☎ 0800 808 855
✉ info@henryschein.co.nz
🌐 www.henryschein.co.nz

HENRY SCHEIN®
DENTAL

Rely on Us™