



RODIN[®]
3D Resin

A Pac-Dent Brand

Glaze N₂-Free

Light-Curable Naturalizing Resin

Glaze N₂-Free is a biocompatible, Class II, light-curable resin designed for surface characterization of dental restorations, including hybrid ceramics, PMMA resins, and CAD/CAM composites. It neutralizes the oxygen inhibition layer for complete polymerization, preventing tacky surfaces and ensuring a smooth, glossy finish that mimics natural enamel.




BENEFITS & FEATURES

- Superior Esthetics
- Sealing and Repairing
- Color Stability
- Versatility
- Intraoral Application



Glaze N₂-Free vs. Optiglaze



*Data provided by Dr. Giordano of Boston University

ORDER INFORMATION:

24031 Glaze N₂-Free, 5 mL Bottle

24015 Glaze N₂-Free, 15 mL Bottle



Pac-Dent, Inc.
670 Endeavor Circle, Brea, CA 92821
www.rodin-3d.com | www.pac-dent.com

Customer Support
(909) 839-0888
info@pac-dent.com

Technical Support
(714) 990-5488
rodin@pac-dent.com



Made in USA