

ahh...way better
GISHY GOO!TM


COMFORT
FOR
EVERY
SMILE

BRACKET
RELIEF
FROM
ULTRADENT



Gishy Goo™ Material

It's the **improved**
replacement for wax.



- Patients say it is easier to place and remove and stays in place longer than wax
- Delicious Bubble Gum flavor

Gishy Goo material is a specially formulated silicone elastomer that can be placed over wires, brackets, fixed appliances, etc. Patients say this unique material is simple to place, stays on longer, and is easier to remove than wax.

Simply squeeze a pea-size amount from the syringe, mix by hand for 10 seconds, and place over the appliance.* There is no messy cleanup and Gishy Goo material will stay in place until it is time to be removed, up to 12 hours!



ULT-500-059 - Gishy Goo **Syringe 10pk**

ULT-2749 - Gishy Goo **Syringe 50pk**
5 ml syringe

* See the IFU for complete instructions

© 2026 Ultradent Products, Inc. All rights reserved.
98999-PKAR01 021626

800.808.855 | ULTRADENT.COM

U2
ULTRADENT