

## 3M™ Impregum™ Super Quick

**Polyether Impression Material** 



# Capture every detail in 2 minutes? Yes, it can.

Now also available in cartridge delivery!



#### Quicker set.

With just a two minute setting time, Impregum Super Quick Material lives up to its name and allows for faster impressions in cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



#### 100% Polyether. Trusted precision.

Impregum Super Quick Materials demonstrate excellent performance in moist situations and exhibit great flow behavior, leading to higher levels of accuracy and detailed, high-precision restorations.

The nature of polyether allows you to get the full 45 seconds of working time, avoiding inaccuracies associated with premature setting times of VPS materials.



# Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth increases your patients' comfort.

### 3M™ Impregum™ Super Quick Polyether Impression Material

- For monophase and 1-step impression techniques.
- Perfect for impressions for singleunit crowns, implants, inlays/onlays and smaller bridges up to 3 units.
- Available in heavy, medium, light and ultra light body consistencies.
- Impregum Super Quick Material is scannable with the current generation of scanners for a more efficient workflow.







## The fastest way to polyether precision for smaller cases.

	Order code	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)	69385	1	Monophase or 1-step	/	~	1	~	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)	69417	7	Monophase or 1-step	~	•	~	<b>✓</b>	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)	69406	1	1-step	<b>*</b>	<b>/</b>	<b>/</b>	<b>✓</b>	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)	69379 (4)	7	1-step	-	•	1	•	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)	69418 (4) 69420 (2)		1-step	-	•	-	•	0:45	2:00

	Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit
69414	3M™ Impregum™ Penta™ Super Quick Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer syringe, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)
69415	3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body  1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 13M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 13M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)
69407	3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body  1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

Refills/Trial Kits
3M™ Impregum™ Penta™ Super Quick Medium Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
3M" Impregum" Super Quick Medium Body Refill 4 cartridges – 50 ml each, 10 3M" Garant™ Mixing Tips Purple
3M <sup>™</sup> Impregum <sup>™</sup> Penta <sup>™</sup> Super Quick Heavy Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
3M™ Impregum™ Super Quick Light Body Refill 4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips Purple
3M™ Impregum™ Super Quick Ultra Light Body Refill 4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips Purple
3M™ Impregum™ Super Quick Ultra Light Body Trial Kit 2 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips Purple

	Trays
71622	3M <sup>™</sup> Dual Arch Impression Trays Posterior, 48 Trays
71623	3M <sup>™</sup> Dual Arch Impression Trays Anterior, 35 Trays
71615	3M <sup>™</sup> Directed Flow Impression Tray, Size S, upper jaw, 10 Trays
71616	3M <sup>™</sup> Directed Flow Impression Tray, Size M, upper jaw, 10 Trays
71617	3M <sup>™</sup> Directed Flow Impression Tray, Size L, upper jaw, 10 Trays
71618	3M <sup>™</sup> Directed Flow Impression Tray, Size S, lower jaw, 10 Trays
71619	3M <sup>™</sup> Directed Flow Impression Tray, Size M, lower jaw, 10 Trays
71620	3M™ Directed Flow Impression Tray, Size L, lower jaw, 10 Trays

	Accessories
69387	3M <sup>™</sup> Impregum <sup>™</sup> Penta <sup>™</sup> Super Quick Medium Body Cartridge
69401	3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge
77580	3M™ Garant™Cartridge Dispenser 1:1/2:1
30600	3M <sup>™</sup> Polyether Tray Adhesive, 1 bottle – 17 ml
69408	3M™ Polyether Contact Tray Adhesive, 1 bottle – 17 ml
77949	3M™ Penta™ Mixing Tips – Red, 50 Tips
71210	3M™ Elastomer Syringe Set, 1 Syringe, 4 Tips, 4 Barrels, 1 Brush
71507	3M™ Intra-Oral Syringes – Purple, 20 Syringes
71508	3M™ Intra-Oral Syringes – Purple, 50 Syringes
71451	3M <sup>™</sup> Polyether Intra-Oral Tips – White, 50 Tips
71454	3M™ Polyether Mixing Tips - Purple, 50 Tips



#### 3M Oral Care

3M Australia Pty Limited Building A, 1 Rivett Road North Ryde, NSW 2113 1300 363 454 3M.com.au/ImpregumSQ **3M New Zealand Limited** 94 Apollo Drive Rosedale, Auckland 0632 0800 808 182 3M.co.nz/ImpregumSQ Now also available in Ultra Light Body.

