

MICROTOUCH® DENTA-GLOVE® Blue Nitrile

DENTA-GLOVE® SERIES EXAMINATION GLOVES

Nitrile, powder-free examination glove with ergonomic design for extra comfort

This non-sterile, nitrile examination glove delivers the ultimate fit and comfort. Designed with ERGOFORM technology, this glove helps reduce muscle effort that can cause stress on joints, ligaments and tendons. Durable enough for rugged donning while providing great tactile sensitivity and grip for fine, detailed procedures.

KEY FEATURES & BENEFITS

- Innovative film technology provides comfort and fit
- Robust strength helps avoid rips and tears while donning and inner polymer coating facilitates ease of donning
- Textured finger surfaces provide a secure grip on instruments in wet or dry conditions
- Tested Safe for Oral Examination, ensuring that the glove is safe for patient contact

RECOMMENDED FOR

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Oral examinations
- Protection from latex Type I allergy in HCW's or patients
- Standard, moderate risk examination procedures

FEATURED TECHNOLOGIES



MICROTOUCH® Everyday, Everywhere Protection

MICROTOUCH® premium medical examination gloves are scientifically developed to fully equip and comfortably protect medical professionals so that they have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICROTOUCH[®] DENTA-GLOVE[®] Blue Nitrile

DENTA-GLOVE[®] SERIES EXAMINATION GLOVES

PRODUCT INFORMATION

Material	Nitrile	Freedom From Holes	1.5 AQL
Not Made From Natural Rubber Latex	Yes	Palm Thickness mm (mils)	0.08 (3.1)
Color	Blue	Allergy Prevention	Latex (Type I)
Cuff Length	Standard	Tested For Use With Chemotherapy Drugs	No
Powder Content	Powder-free		
External Glove Surface	Textured Fingers	Sterile	No

ORDERING INFORMATION

Size	X-Small	Small	Medium	Large	X-Large
Ansell Product Code	3130050XS	31300500S	31300500M	31300500L	3130050XL

PACKAGING AND STORAGE

Packaging	100 gloves per box (size XS-L); 90 gloves per box (size XL); 10 boxes per carton/case
Shelf Life	3 years
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Medical Device: Class I, EN 455 Part 1, EN 455 Part 2, EN 455 Part 3, EN 455 Part 4, ISO 13485 :2003, ISO 9001 :2008, ISO 10993-10 :2002



NO RUBBER LATEX

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Latin America & Caribbean
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
 Brazil Tel: +55-11-3356-3100
 Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0)2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com