

MICRO-TOUCH® AFFINITY™

Performance Series

Neoprene, powder-free examination glove with excellent elasticity and barrier protection

This non-sterile neoprene synthetic formulation provides superior protection, comfort, fit, elasticity and tactile sensitivity with excellent resistance to a broad range of chemicals.

KEY FEATURES AND BENEFITS:

- · Offers comfort level similar to latex.
- Textured finger surfaces provide a secure grip in both dry and wet conditions for enhanced instrument handling.
- Excellent tear and puncture resistance
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs*
- * See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Enhanced grip in wet environments
- Sustained, high-risk examination procedures

FEATURED TECHNOLOGIES







MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



MICRO-TOUCH® AFFINITY™

Material	Neoprene (Polychloroprene)		Not Made From Natural Rubber Latex		Yes	Yes	
Color	Green			f Length	Standard	Standard	
Powder Content	Powder-Free			ernal Glove Surface	Textured Fi	Textured Fingers	
Freedom from Holes (Inspection level I)	1.5 AQL		Palm Thickness (mm/mil)		oil) 0.13 / 5.1	0.13 / 5.1	
Finger Thickness (mm/mil)	0.16 / 6.3			ergy Prevention	Latex (Type	Latex (Type I)	
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)		Tested For Use With Chemotherapy Drugs			Yes, in accordance with ASTM D6978 and US FDA cleared	
Sterile	No			istatic	Not Tested	Not Tested	
Glove Length (mm/inches)	240 / 9.5		Product Segmentation		Performanc	Performance Series	
ORDERING INFORMATIO	N						
Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	
Product Code	3770	3771		3772	3773	3774	
PACKAGING AND STORA	GE						

STANDARDS AND CERTIFICATIONS

5 Years

ignition.

Packaging Shelf Life

Storage Instructions

ASTM D6977, CA Prop 65, Category III, EN 455 Part 1, EN 455 Part 3, FDA21 CFR 177-2600, ISO 13485, ISO 10993-10:2002

100 gloves per dispenser 10 dispensers per case 1000 gloves per case

Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 1-844-494-7854 insidesalesus@ansell.com Latin America & Caribbean Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa Tel: +6 Tel: +32 (0) 2 528 74 00 Email: Email: info@ansell.eu

Asia Pacific Tel: +603 8310 6688 apac.medical@ansell.com Australia & New Zealand Tel: +61 3 9270 7270 protection@ap.ansell.com