

MICRO-TOUCH® DENTA-GLOVE® NRL

DENTA-GLOVE® Series

Latex, powder-free, examination glove with maximum comfort for a wide range of dental procedures

This non-sterile glove is an ideal choice for dental specialists who switch gloves frequently and experience hand fatigue from continuous use.

KEY FEATURES AND BENEFITS:

- Latex formulation and optimized film thickness offer a high level of elasticity, flexibility. and sensitivity
- · Inner chlorinated surface to ensure fast and easy donning
- Tested Safe for Oral Examination, ensuring that the glove is safe for patient contact

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

Recommended For

- Enhanced grip in wet environments
- Standard, moderate risk examination procedures
- Oral examinations





MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



MICRO-TOUCH® DENTA-GLOVE® NRL

PRODUCT INFORMATION				
Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No	
Color	Natural, White	Cuff Length	Standard	
Powder Content	Powder-Free	External Glove Surface	Fully Textured	
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.14 / 5.5	
Finger Thickness (mm/mil)	0.16 / 6.3	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	
Tested For Use With Chemotherapy Drugs	No	Sterile	No	
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5	
Product Segmentation	DENTA-GLOVE® Series			

ORDERING INFORMATION					
Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Ansell Product Code	04653	04654	04655	04656	04657

PACKAGING AND STORAGE			
Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *90 gloves per dispenser (Size XL) *900 gloves per case (Size XL)		
Shelf Life	3 Years		
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.		

STANDARDS AND CERTIFICATIONS

ASTM D3578, Category III, EN 420:2003 + A1:2009, EN 455 Part 1, EN 455 Part 2, EN 455 Part 3, EN 455 Part 4, ISO 11193-1:2008, ISO 9001, ISO 13485, ISO 10993-10:2010

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 1-844-494-7854 insidesalesus@ansell.com Latin America & Caribbean Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa Tel: +6 Tel: +32 (0) 2 528 74 00 Email: Email: info@ansell.eu

Asia Pacific Tel: +603 8310 6688 apac.medical@ansell.com Australia & New Zealand Tel: +61 3 9270 7270 protection@ap.ansell.com