## TISSUE MANAGEMENT

## Astringedent™

15.5% FERRIC SULFATE









- The "Classic" hemostatic agent for profound hemostasis
- Stops moderate bleeding
- Eliminates sulcular fluid contamination for optimal bonding
- · Decreases costly impression remakes

Astringedent hemostatic is an aqueous 15.5% ferric sulfate hemostatic solution with a pH of  $\sim$ 1.0.

Astringedent hemostatic solution is well suited for a variety of dental and oral surgery procedures to arrest bleeding. Astringedent hemostatic can be used to prevent leakage caused by sulcular fluid contamination during direct bonding procedures.

Listed as a "CAN'T LIVE WITHOUT" product by a prominent independent research institute.<sup>2</sup>

Note: ViscoStat<sup>™</sup> and Astringedent hemostatic agents should be used with a Metal Dento-Infusor<sup>™</sup> tip. The plastic Blue Mini<sup>™</sup> Dento-Infusor<sup>™</sup> tip should be used when you are dealing with newly healed epithelium, as the softer tip is slightly less aggressive.



0111 - Astringedent Bottle 1pk 0686 - Astringedent IndiSpense™ Syringe 1pk

1. realityesthetics.com. 2. "Can't Live Without" Clinical Research Associates Newsletter, Volume 21, Issue 7, July 1997.

## Astringedent™ X

12.7% IRON SOLUTION







- Clinicians "go-to" hemostatic for all case situations
- Stops minor to severe bleeding

Astringedent X hemostatic is an aqueous 12.7% iron solution that works quickly to stop bleeding in seconds. It contains equivalent ferric sulfate and ferric subsulfate. Note: Diluted Astringedent X hemostatic does not equal ViscoStat or Astringedent hemostatics.

Use when a stronger, more potent hemostatic is required and when the attainment of quality hemostasis may be more challenging (e.g., in cases of difficult-to-stop, problem bleeding).





Astringedent X hemostatic and Metal Dento-Infusor tip facilitate profound hemostasis, even with challenging cases.





0112 - Astringedent X Bottle *1pk* 690 - Astringedent X <u>IndiSpense Syringe</u> *1pk* 30 ml

## **Astringedent™ Spot Remover**

**CLEANING SOLUTION** 

Astringedent Spot Remover is designed to remove ViscoStat hemostatic, Astringedent hemostatic, and Astringedent X hemostatic stains that will not come out of clothing with soap and water. Not for intraoral use.

2160 - Astringedent Spot Remover 1pk 30 ml bottle

