



Futura

Designed to make life easier, reduce costs and guarantee higher efficiency

- Can be used with Tap Water. Utilises an in-Built Demin Filter so that Bottled Demin water or external filtration systems are a thing of the past. Saves on water costs and reduces environmental impact
- 4.3 inch colour touch display. For simplified function navigation.
- LED light bar to illuminate the workspace in front of the Steriliser.
- Heavy Duty vacuum pump for speed and reliability
- Optional Barcode and cycle report printer. Track critical instruments digitally.
- Delayed start so Vacuum and Helix tests can be programmed to run before a practice opens. A real time saver.
- Wi-Fi & Ethernet port are standard inclusions so Steriliser can be networked.
- User management system allows users to be tracked to cycle data via personal PIN codes. Complies with Australian and NZ standards and means no paper printout is required.
- Cycle data saved as a PDF file meaning no software required to save data. Can be downloaded to PC via included USB.
- Steriliser can be installed as a stand alone device or connected to mains water or external R/O systems. Multiple options to cover varying installation requirements.
- Built in water conductivity sensor to ensure water quality is always at it's best.
- Optional Cloud based software to manage machine performance, cycle data and provide remote access for technical departments.

Industry leading 3 year or 3,000 cycle warranty!