

Bioclear® Daily DISINFECTANT

EDR DENTAL LINIT WATER LINES



Bioclear Daily is a disinfectant liquid for continuous use in dental unit water lines (DUWLs).

Used daily, it reduces the planktonic bacteria count in water flowing through the unit, minimising the potential for biofilm formation.

Pack Size:

1L Bottle, including measuring cap for accurate dosage.

Key Points:

- For continuous use in dental unit waterline feed bottles.
- Continuously controls water quality to reduce the potential formation of biofilm.
- Economical as it is used at 1% concentration ie only 10ml per 1L.
- Prepared solutions can be kept insuitable plastic containers for up to 10 days.
- Excellent material compatibility.
- Simple and easy to use.
- Can also be used in dental units with compatible built-in automated dosing systems. (Ensure the dosing equipment is set to deliver a concentration of 1%.)

Item code: QD-DS0241





