

Wireless Endo Motor with Apexlocator

Powerful functions,

Precise measurements,

Efficient and safe,



Endo motors



Order name: C-SMART miniAP

This product has been applied layout-design patent, Counterfeit must be investigated.





Feature:

- 1-Motor and Apex locator two in one
- 2-Integrated mini head contra angle
- 3-Fine color display screen-Building muilt file system
- 4-Auto self-test and calibration-Auto switch work model
- 5-Adjustable angle and adjustable speed of reciprocating motion
- 6-Contact-less charging
- 7-Large capacity battery
- 8-Switchable display for right hand or left hand use
- 9-The root electrode is built in the head, and the root canal measure ment function can be well realized without the external root mea surement buckle.

Technical:

Adaptor working voltage	100-240Vac,50/60HZ;
Adjustable torque,range	0.6-3.9N.cm(6-39mN.m)
5 setting for adjustable rotational speed ranging:	125rpm-600rpm
Volume(cm)	19X13X6
Weight(kg)/PC\$	0.7
Packing material	carton
QTY per master carton	12 pc/ctn







• auto reverse by torque • auto reverse by apical position



