



Performance Meets Comfort.  
Designed for Patients,  
Built for Surgeons.

CyberTech  
**SURGICAL  
SUTURES**

## CYBERTECH SURGICAL SUTURES

### FEATURES

PGA (Polyglycolic Acid) absorbable suture is a synthetic polymer that naturally dissolves within weeks\*






- Predictable absorption for efficient wound support
- Secure needle-to-suture attachment engineered for strength, consistency, and reduced tissue trauma
- Long cutting-edge design for precise handling and suturing control
- Single-use, sterile



\*Refer to the instructions for use for detailed information.






### PGA FA – Ideal for short-term wound support

Designed for cases requiring fast absorption, PGA FA provides reliable, temporary wound closure during the early healing phase.

ITEM CODE	DESCRIPTION	TYPE	COLOR	USP <sup>1</sup>	NEEDLE CURVATURE	NEEDLE LENGTH	THREAD LENGTH	NEEDLE TIP	NEEDLE SHAPE	POINT GEOMETRY	THREAD
9008112	Synthetic, Coated (Not Antimicrobial)	Absorbable (Fast)	Undyed	3/0	$\frac{3}{8}$ Circle	18mm	45cm	Reverse Cutting	Triangle		Braided
9008113	Synthetic, Coated (Not Antimicrobial)	Absorbable (Fast)	Undyed	4/0	$\frac{3}{8}$ Circle	18mm	45cm	Reverse Cutting	Triangle		Braided
5723794	Synthetic, Coated (Not Antimicrobial)	Absorbable (Fast)	Undyed	4/0	$\frac{1}{2}$ Circle	19mm	45cm	Taper Point	Round		Braided
5723795	Synthetic, Coated (Not Antimicrobial)	Absorbable (Fast)	Undyed	3/0	$\frac{1}{2}$ Circle	19mm	45cm	Taper Point	Round		Braided
5723796	Synthetic, Coated (Not Antimicrobial)	Absorbable (Fast)	Undyed	3/0	$\frac{1}{2}$ Circle	22mm	45cm	Taper Point	Round		Braided

### PGA – Ideal for longer-term healing needs

Offers extended tensile strength and slower absorption, making it suitable for wounds that require prolonged support.

ITEM CODE	DESCRIPTION	TYPE	COLOR	USP <sup>1</sup>	NEEDLE CURVATURE	NEEDLE LENGTH	THREAD LENGTH	NEEDLE TIP	NEEDLE SHAPE	POINT GEOMETRY	THREAD
9008114	Synthetic, Coated (Not Antimicrobial)	Absorbable	Violet	4/0	$\frac{1}{2}$ Circle	18mm	45cm	Taper Point	Round		Braided
9008115	Synthetic, Coated (Not Antimicrobial)	Absorbable	Violet	4/0	$\frac{3}{8}$ Circle	18mm	70cm	Reverse Cutting	Triangle		Braided
9008116	Synthetic, Coated (Not Antimicrobial)	Absorbable	Violet	3/0	$\frac{3}{8}$ Circle	24mm	70cm	Reverse Cutting	Triangle		Braided
5723797	Synthetic, Coated (Not Antimicrobial)	Absorbable	Violet	3/0	$\frac{1}{2}$ Circle	19mm	45cm	Taper Point	Round		Braided
5723798	Synthetic, Coated (Not Antimicrobial)	Absorbable	Violet	3/0	$\frac{1}{2}$ Circle	22mm	70cm	Taper Point	Round		Braided

<sup>1</sup>USP: United States Pharmacopeia

Order now from [HenrySchein.co.nz](https://henryschein.co.nz)