



# EDGEAPEX HP™

**EDGEAPEX HP™** Designed to enhance precision and efficiency in endodontic procedures. Cordless and with an Integrated Apex Locator for maximum convenience.

## Benefits

- Compatible with the majority of the existing file systems for better adaptability.
- Exceptional versatility with two movement options, continuous rotation and reciprocation.
- Cordless portable motor with wireless charging, auto on/off, and an integrated Apex Locator for maximum convenience.
- High-performance brushless motor with quiet operation and an extended service life.
- It adopts real-time feedback technology and dynamic torque contra angle, effectively preventing file separation.

## Specifications:

- Great flexibility with a 360° rotation contra angle.
- Ultra-small head and slim neck for enhanced visibility and improved access to difficult areas.
- 10 available memories and pre-program settings for EdgeEndo file systems to save time.
- 4 available motions: CW / CCW / REC / ATR
- Multiple Apical functions
- Right-handed and left-handed adjustable display
- RPM range: 100 - 2500
- Torque range: 0.4 - 5.0 Ncm



## The box includes

A motor handpiece, a charging base, a contra angle, a measuring wire, lip hook (x2), file clip (x4), touch probe (x2), power adapter, protective silicon cover (x2), spray nozzle, o-rings (x2), disposable sleeves (1 pack).

Motor Hand piece: **EN-EAPEXRU**  
Charging Base: **EN-EAPEXCB**

Item Code: **EN-EAPEXHP**

## Main accessories

Product	Item Code	Image
<b>MEASURING WIRE</b>	EN- EAPEXNMW	
<b>TOUCH PROBE</b>	EN- EAPEXTP	
<b>LIP HOOK</b>	EN- EAPEXLH	
<b>CONTRA ANGLE HANDPIECE CA 161</b>	EN- EAPEXCA	
<b>FILE CLIP</b>	EN- EAPEXFC	

**Warranty:** 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.