



EDGEAPEX™

EDGEAPEX™ Effectively locates the apical constriction, enabling precise determination of root canal working length for optimal treatment outcomes.



Main Unit: **EN-EEALURU**
Charging Base: **EN-EEALUB**

Benefits

The EdgeApex is an apex locator and a pulp vitality tester.

It utilizes cutting-edge multi-frequency technology combined with a highly efficient digital signal processor, delivering optimal accuracy and featuring a high resolution large screen.







Pulp vitality tests diagnose pulpitis, determine the need for root canal treatment or other dental procedures, and assess overall tooth health.

The box includes

- Charging Base (1)
- Measuring Wire (1)
- File Clip (4)
- Lip Hook (5)
- Touch Probe (2)
- Pulp Tester Probe (2)
- Tester (1)
- Lithium-ion Battery (1)
- Power Adapter (1)
- Enlarged display of apical foramen (1)

Item Code: **EN-EEALU**

Main Accessories

Product	Item Code	Image
MEASURING WIRE	EN-EEALUMW	
TOUCH PROBE	EN-EEALUMF	
LIPHOOK	EN-EEALULH5	
ENDO RULER	EN-EEALUER	
FILE CLIP	EN-EEALUFC5	
PULP TESTER PROBE	EN-EEALUPTP	

Warranty: 1 year for device. 6 months for accessories.