



# EDGEFLOW™

360° Rotating with excellent gutta-percha fluidity for a uniform, and smooth, bubble-free output ensuring seamless fillings.

## Benefits

- Excellent gutta-percha fluidity for a uniform, bubble-free output ensuring seamless fillings.
- Wireless and electric design for maximum comfort and ease of use.
- Dual battery system provides extended endurance, allowing for approximately 100 root canal fillings.
- Precise and intuitive display to better monitor the remaining gutta-percha.

### Specifications:

- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Reusable injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast



Main Unit: **EN-EEPTORU**  
Charging Base: **EN-EEPTOCH**

## The box includes

A main unit, Charging Stand, Power Adapter, Battery (x2), Heat Insulation Hood, Wrench of Needle, Gutta-Percha Ring, Injection Needle (x4) E23G-M (1) - E23G-L (1) - E25G-L (2).

Item Code: **EN-EEPTO**

## Main Accessories

Product	Item Code	Image
<b>WRENCH</b>	EN-EEPTOW	
<b>LITHIUM-ION BATTERY</b>	EN-EEPTOB	
<b>THERMAL PROTECTOR CAP</b>	EN-EEPTOHIH3	

## Tips

Product	Item Code	Image
<b>E20G-S INJECTION NEEDLE</b>	EN-EE20GSIN	
<b>E20G-M INJECTION NEEDLE</b>	EN-EE20GMIN	
<b>E20G-L INJECTION NEEDLE</b>	EN-EE20GLIN	
<b>E23G-M INJECTION NEEDLE</b>	EN-EE23GMIN	
<b>E23G-L INJECTION NEEDLE</b>	EN-EE23GLIN	
<b>E25G-M INJECTION NEEDLE</b>	EN-EE25GMIN	
<b>E25G-L INJECTION NEEDLE</b>	EN-EE25GLIN	

**Warranty:** 1 year for device.