



EDGEPACK™

EDGEPACK™ Features high-precision temperature control with pre-set working temperatures, ensuring optimal compatibility with various gutta-percha types for precise and reliable obturation.

Benefits



- Wireless Gutta-Percha Obturation Device for seamless, cordless operation.
- Ergonomic and ambidextrous design for enhanced comfort and usability.
- 0.35mm pre-bend tip for improved access and better penetration into narrow root canals.
- High-precision temperature control for consistent and reliable obturation.
- Efficient GP cutting system.
- Compact size for easy cleaning and convenient storage.
- Visual and auditory prompts for enhanced safety during procedures.

Specifications:

- Preset temperature options: 150°C, 180°C, 200°C, and 230°C.
- Advanced time-related features:
 - Ultra-fast warm-up time of 0.2 seconds.
 - Auto shut-off after 10 seconds for energy efficiency.
 - Long-lasting battery life—provides up to 4 hours per charge.


The box includes

A Charging Stand, Power Adapter, WP4004 Tip, WP4504 Tip, WP5506 Tip, WP5508 Tip, Heat Insulation Sleeve of Tips (x4).


Main Unit: **EN-EEDPURU**





Item Code: **EN-EEDPU**

Main Accessories

Product	Item Code	Image
CHARGING BASE	EN-EEDPUCS	

Tips

Product	Item Code	Image
WP4004 WORK TIP	EN-EEDPUT4004	
WP4504 WORK TIP	EN-EEDPUT4504	
WP5004 WORK TIP	EN-EEDPUT5004	
WP5506 WORK TIP	EN-EEDPUT5506	

Product	Item Code	Image
WP5508 WORK TIP	EN-EEDPUT5508	
WP5508L WORK TIP	EN-EEDPUT5508L	
WP5510 WORK TIP	EN-EEDPUT5510	
WP6012 WORK TIP	EN-PUT6012	

Warranty: 1 year for device and tips.