Fit Strip Flexible

I nterproximal finishers with a

Twist



FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts,

Garrison Dental Solutions.

One handy handle! FitStrips' unique attachable handle has been IMPROVED! It's even easier to attach (or remove) and makes a world of difference for both clinician ease-of-use and patient comfort. Your hand stays outside the mouth, providing a clear field of view and easier operation.

Do the Twist! Twist the colourcoded barrel to adjust FitStrips' curvature and it's automatically locked in place—curved for finishing or straight for IPR. It doesn't get any easier than that!

FitStrip™ – make your procedures easier!

- Orthodontics: The first truly complete IPR/tooth slenderising kit
- Composite restorations: contouring/finishing all interproximal surfaces
- Crown and bridge: cement removal and clean-up

Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

