Fit Strip.

I nterproximal finishers with a



FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, Garrison Dental Solutions.

One handy handle! FitStrips' unique attachable handle has been IMPROVED! It's even easier to attach (or remove) and makes a world of difference for both clinician ease-ofuse and patient comfort. Your hand stays outside the mouth, providing a clear field of view and easier operation.

Do the Twist! Twist the colourcoded barrel to adjust FitStrips' curvature and it's automatically locked in place—curved for finishing or straight for IPR. It doesn't get any easier than that!

FitStrip™ - make your procedures easier!

- Orthodontics: The first truly complete IPR/tooth slenderising kit
- Composite restorations: contouring/finishing all interproximal surfaces
- Crown and bridge: cement removal and clean-up

Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.



Fit Strip.



FitStrip™

Flexible

Interproximal finishers with a

Twist

Your Complete IPR Solution!

NEW! .40mm and .50mm strips easily accommodate maximum reduction cases

NEW! IPR feeler gauges for exact reduction

IMPROVED! Attachable handle for greater ease of use and patient comfort

Kits All kits come in a package of 10 FitStrips with 2 handles.

NEW Universal Kit

1 IPR gauge

4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 6 double-sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)

10 Strips (2 each: .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

GA1-FPSK04

Refills		Single Sided	Double Sided
Super fine 15 Micron		GA1-FPSXFSS 0.08 mm thick	GA1-FPSXFDS 0.11 mm thick
Fine 30 Micron		GA1-FPSFSS 0.10 mm thick	GA1-FPSFDS 0.15 mm thick
Medium 46 Micron		GA1-FPSMSS 0.13 mm thick	GA1-FPSMDS 0.21 mm thick
Coarse 90 Micron		GA1-FPSCSS 0.18 mm thick	GA1-FPSCDS 0.30 mm thick
Medium 46 Micron			GA1-FPSTMDS 0.40 mm thick
Medium 46 Micron			GA1-FPSXTMDS 0.50 mm thick
Serrated		GA1-FPSSAW 0.05 mm thick	
FitStrip replacement handles (2-pack) GA1-FPSHANDLE-2			
Fitstrip IPR Gauge GA1-FPSGAUGE			A1-FPSGAUGE



