



FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts,

Garrison Dental Solutions.

One handy handle! FitStrip's unique attachable handle has been IMPROVED! It's even easier to attach (or remove) and makes a world of difference for both clinician ease-of-use and patient comfort. Your hand stays outside the mouth, providing a clear field of view and easier operation.

Do the Twist! Twist the colorcoded barrel to adjust FitStrips' curvature and it's automatically locked in place—curved for finishing or straight for IPR. It doesn't get any easier than that! See our video online!

## FitStrip™ makes your procedures easier!

- Othodontics: The first truly complete IPR/ tooth slenderizing kit
- Composite restorations: Contouring/ finishing all interproximal surfaces
- Crown and bridge: Cement removal and clean-up

Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.





Flexible
I nterproximal finishers with a
Twist



**NEW!** .40mm and .50mm strips easily accommodate maximum reduction cases

**NEW!** IPR gauges for exact reduction


**IMPROVED!** Attachable handle for greater ease of use and patient comfort

Kits All kits come in a package of 10 FitStrips with 2 handles

This fill kits come in a package of to this tips with 2 handles.		
NEW Universal Kit for IPR and finishing/polishing		
1 IPR gauge	CA1 FDCKO4	
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	GA1-FPSK04	
6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)		
Starter Kit for finishing/polishing		
2 serrated saw	CA1 FDCKO1	
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	GA1-FPSK01	
4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)		

Refills	Single Sided	Double Sided
Super fine 15 Micron	GA1-FPSXFSS 0.08 mm thick	GA1-FPSXFDS 0.11 mm thick
Fine 30 Micron	GA1-FPSFSS 0.10 mm thick	GA1-FPSFDS 0.15 mm thick
Medium 46 Micron	GA1-FPSMSS 0.13 mm thick	GA1-FPSMDS 0.21 mm thick
Coarse 90 Micron	GA1-FPSCSS 0.18 mm thick	GA1-FPSCDS 0.30 mm thick
Medium 46 Micron		GA1-FPSTMDS 0.40 mm thick
Medium 46 Micron		GA1-FPSXTMDS 0.50 mm thick
Serrated	GA1-FPSSAW 0.05 mm thick	
FitStrip™ replacement handles (2-pack) GA1-FPSHANDLE-		GA1-FPSHANDLE-2
FitStrip™ IPR Gauge		GA1-FPSGAUGE







Prices exclude GST and freight (E&OE).

Pricing not applicable to Govt or Corporate Accounts.

To order telephone 0800 808 855, 8am to 5pm | freefax 0508 808 555, 24 hours orders @henryschein.co.nz | www.henryschein.co.nz