

When design and function come together...you get **FUSION!**

Composi-Tight® 3DFusion™ Sectional Matrix System



Tight contacts – guaranteed!

Superior tooth separation force through advanced engineering and material selection guarantees perfect contacts every time.

Ideal anatomy!

Soft-Face™ silicone ring tips with enhanced marginal ridge anatomy hug the pre-contoured dead-soft matrices to the tooth for flash free, perfectly contoured restorations.

Simplified technique!

Composi-Tight® 3D Fusion™ is packed with features like easy ring stacking, over-the-wedge placement and Ultra-Grip™ anti-ring slip extensions that make system use intuitive and frustration free.

Matrix band perfection



Composi-Tight® 3D Fusion™ Full Curve Non-Stick Matrices are three dimensionally contoured, ultra-thin and dead-soft to help you recreate ideal interproximal anatomy.



short ring



tall ring




wide prep ring

It just keeps going...

Enhanced durability for extended ring life and greater tooth separation

Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



Fins with serious function

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.

Wider indications, wider smiles!




Dentistry's first Wide Preparation separator ring (green) makes Composi-Tight® 3D Fusion™ the most versatile and user-friendly sectional matrix system available. The Wide

Prep ring greatly simplifies what until now has been one of the most challenging posterior composite restorations. Now you can restore these big preps with confidence knowing that you will achieve excellent contact and contour without having to take extra steps.

Only a sectional matrix can deliver results like this!

Tofflemire-style System

- ❌ Fails to restore proximal anatomy
- ❌ Thin contact at the marginal ridge
- ❌ Large food trap below
- ❌ Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Composi-Tight® 3D Fusion™ Sectional Matrix System

- ✓ Operator-friendly retaining system
- ✓ Naturally contoured bands
- ✓ Anatomically correct contacts
- ✓ Contacts at the height of contour
- ✓ Rings produce optimal tooth separation for excellent, tight contacts



Composi-Tight® 3DFusion™

Sectional Matrix System



Kit FX-MFF-00

Contents:

- 1 – Short Ring (blue)
- 1 – Tall Ring (orange)
- 25 – Assorted Matrix Bands
- 20 – Assorted Ultra-Adaptive Wedges
- 1 – Ring Placement Forceps



Kit FX-KFF-50

Contents:

- 1 – Short Ring (blue)
- 1 – Wide Prep Ring (green)
- 40 – Assorted Matrix Bands
- 40 – Assorted Ultra-Adaptive Wedges
- 1 – Ring Placement Forceps

Includes dentistry's premier heavy duty stainless steel ring placement forceps, which place virtually any ring.

For more information go to www.garrisdental.com. Be sure to check out the website to see the full line of innovative Garrison products.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456



Mention Code
INNZHS19