


Composi-Tight® 3D Fusion™

Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.

Get a grip The universal gripping block gives you absolute control of the wedging process whether you're using a cotton pliers, hemostat or "pin" tweezers. Use what you like – Composi-Tight® 3D Fusion™ has you covered.



A real softie Composi-Tight® 3D Fusion's Soft-Face™ over-mold allows the wedge to do what no other wedge can truly do – actually adapt to interproximal irregularities.

Tough on the inside

The firm inner skeleton allows for easy insertion and tooth separation just like a traditional wedge.

No pain – all gain

The upturned ski-like tip and soft radius edges glide over the rubber dam or papilla for trauma free insertion.




Fins with serious function

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.

Keep it clean (and dry!)

Composi-Tight® 3D Fusion's patented design easily follows the contours of the teeth sealing things up nice and tight while preventing overhangs.



Description	Item	Qty.
Composi-Tight® 3D Fusion™ Wedge Kit	FXK4	400
Composi-Tight® 3D Fusion™ Wedge Kit	FXK4-M	200

Refills	Standard		Bulk	
	Item	Qty.	Item	Qty.
X-Small Yellow	FXYL-M	50	FXYL	100
Small Blue	FXBL-M	50	FXBL	100
Medium Orange	FXOR-M	50	FXOR	100
Large Green	FXGR-M	50	FXGR	100



Composi-Tight® 3D Fusion™ Wedges fuse two materials together to produce the perfect combination of adaptability and retention.

Garrison
Dental Solutions

Composi-Tight® 3D Fusion™ Ultra Adaptive Wedges




Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.

Composi-Tight® 3D Fusion™ Wedge Kits

Description	Item	Qty.
Composi-Tight® 3D Fusion™ Wedge Kit	FXK4	400
Composi-Tight® 3D Fusion™ Wedge Kit	FXK4-M	200

Refills	Item	Standard		Bulk	
		Qty.	Item	Qty.	
X-Small Yellow	FXYL-M	50	FXYL	100	
Small Blue	FXBL-M	50	FXBL	100	
Medium Orange	FXOR-M	50	FXOR	100	
Large Green	FXGR-M	50	FXGR	100	



WedgeWands®

WedgeWands® also available in Cure-Through


WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.

WedgeWands®

Description	Opaque / Clear	Qty.
WedgeWands Kit	WK4 / WCK4	400

Refills	Opaque / Clear	Standard		Bulk	
		Qty.	Opaque / Clear	Qty.	
X-Small Yellow	SWYL / SWCYL	100	WYL / WCYL	300	
Small Blue	SWBL / SWCBL	100	WBL / WCBL	300	
Medium Orange	SWOR / SWCOR	100	WOR / WCOR	300	
Large Green	SWGR / SWCGR	100	WGR / WCGR	300	

A+Wedge® Controls Gingival Bleeding




- Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva
- Improves clarity of restoration area
- Eliminates need to apply an astringent separately

A+Wedge® Kits

Description	Item	Qty.
A+Wedge Kit	GWAK4	200

Refills	Item	Bulk	
		Qty.	Item
X-Small Yellow	GWAYL	100	WYL / WCYL
Small Blue	GWABL	100	WBL / WCBL
Medium Orange	GWAOR	100	WOR / WCOR
Large Green	GWAGR	100	WGR / WCGR



G-Wedge®

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.

G-Wedge™

Description	Item	Qty.
G-Wedge Kit	GWK4	400

Refills	Item	Standard		Bulk	
		Qty.	Item	Qty.	
X-Small Yellow	GSWYL	100	GWYL	300	
Small Blue	GSWBL	100	GWBL	300	
Medium Orange	GSWOR	100	GWOR	300	
Large Green	GSWGR	100	GWGR	300	

SoftWedge®



Soft Wedges™ compress during insertion, then swell to help fill voids to ensure that your preparations stay clean and dry—an absolute must for today's bonded restorations!

Soft Wedge™

Description	Item	Qty.
Soft Wedge Kit	WDK4	400

Refills	Item	Standard		Bulk	
		Qty.	Item	Qty.	
X-Small Yellow	WSDYL	100	WDYL	300	
Small Blue	WSDBL	100	WDBL	300	
Medium Orange	WSDOR	100	WDOR	300	
Large Green	WSDGR	100	WDGR	300	

Be sure to check out the website at www.garrisdental.com for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

Offer expires 12/31/2020
To order: Call toll free **888.437.0032**
or online at www.garrisdental.com



Mention Code
INAU120FX