

Zero in on your target.

Apex ID[®] Digital Apex Locator

Apex ID is compact and modern with unique functions which allow you to focus on what's really important; providing the best care for your patients.

User-set zero point

Customize the "zero point" based on your clinical judgment. No need to constantly monitor the visual display, the audible signal will let you know when you've reached your preferred zero point.

Accuracy in almost any condition

Apex ID has been designed to work in almost any canal condition. Whether dry, wet, bleeding, filled with saline, EDTA, NaOCl, chlorhexidine, the Apex ID can still provide accurate, consistent and dependable readings.

Continuous calibration

Unlike most locators on the market today, Apex ID continuously calibrates for accurate, uninterrupted readings in real-time.




Actual Size
(2.56" x 3.74" or
6.5cm x 9.5cm)

Easy to read
color display

Optional
user-set
zero point

Cordless
power

Weighted snap-in
base for stability



Part #	Description
972-0090	(1) Apex ID, (1) Apex ID Stand, (5) Lip Hooks, (2) File Holder B, (1) File Holder A (1) Probe Cord, (1) Apex ID IFU, (3) AAA 1.5V Batteries
972-0092	(5) Apex ID Lip Hooks
972-0093	(2) Apex ID File Holder B
972-0094	(1) Apex ID File Holder A
972-0095	(1) Apex ID Probe Cord

077-3145 Rev. A