






# Strong, Safe and Versatile.


-  **K-FLEX FILES**
-  **TRIPLE-FLEX FILES**
-  **K-FILES**
-  **HEDSTROM FILES**
-  **REAMERS**

## Stainless Steel Files

Open the path to Endodontic Excellence


### K-Flex Files

K-Flex® Files feature a unique rhomboid design, which provides a short axis for flexibility while maintaining cross-sectional strength. K-Flex is more flexible than conventional "K" files of comparable size, particularly in sizes #30 and above, which facilitates the negotiation of curved canals and their preparation with minimal potential for ledging and perforating.




### Triple-Flex Files

Triple-Flex® Files are triangular-bladed, stainless steel instruments designed for flexibility with integrity. The twisted flutes allow aggressive cutting and debris removal you demand without sacrificing the instrument's inherent stability. Triple-Flex Files are highly breakage resistant, even in the most stressful situations.



### K-Files

The K-File is the strongest of the handfiles with less separation potential. K-Files can bypass obstructions with far greater ease. They are designed to provide a smooth tactile feel during instrumentation.



#### K-Flex Files (pkg of 6)

Size	21mm	25mm	30mm
80		15174	15258
70		15170	15254
60	15334	15166	15250
55	15330	15162	15246
50	15326	15158	15242
45	15322	15154	15238
40	15318	15150	15234
35	15314	15146	15230
30	15310	15142	15226
25	15306	15138	15222
20	15302	15134	15218
15	15298	15130	15214
10	15294	15126	15210
08	15290	15122	15206
06	60120	60144	60193
15-40 Asst	821-4021	821-4025	821-4030
45-80 Asst	821-4125		821-4130

#### Triple-Flex Files (pkg of 6)

Size	21mm	25mm	30mm
80	25941	25957	25973
70	25940	25956	25972
60	25939	25955	25971
55	25938	25954	25970
50	25937	25953	25969
45	25936	25952	25968
40	25935	25951	25967
35	25933	25949	25965
30	25932	25948	25964
25	25931	25947	25963
20	25930	25946	25962
15	25929	25945	25961
10	25928	25944	25960
08	25927	25943	25959
15-40 Asst	821-3021	821-3025	821-3030
45-80 Asst	821-3121	821-3125	821-3130

#### K-Files (pkg of 6)


Size	21mm	25mm	30mm
80	06026	06099	06170
70	06024	06085	06166
60	06023	06077	06163
55	06022	06076	06162
50	06021	06075	06161
45	06020	06074	06160
40	06019	06069	06159
35	06016	06067	06157
30	06015	06066	06156
25	06014	06065	06155
20	06005	06064	06153
15	06004	06063	06152
10	06003	06062	06151
08	14127	14128	14129
06	62249	62250	62251
15-40 Asst	821-8021	821-8025	821-8030
45-80 Asst	821-8121	821-8125	821-8130

To order call SybronEndo at 800.346.3636, contact your local representative or your authorized SybronEndo dealer.

# Stainless Steel Files

## Pathfinder

Stainless steel Pathfinder™ instruments feature a minimal taper and are used for negotiating difficult or calcified canals.




### Pathfinder (pkg of 6)

21mm	25mm
17023	17017

## Pathfinder CS

Pathfinder® CS instruments feature the same minimal taper as stainless steel Pathfinders but are made with carbon steel to help negotiate difficult or calcified canals faster and more efficiently.




### Pathfinder CS (pkg of 6)

Size	Handle	19mm	21mm	25mm
K1	Dark Brown	22935	22936	23465
K2	Orange	22937	22938	23467

## Hedstrom Files

Hedstrom Files are ground from stainless steel wire. The flute depth increases with the taper, creating a sharper, more flexible instrument.




### Hedstrom Files (pkg of 6)

Size	21mm	25mm	30mm
70	12796	11354	13278
60	12793	11353	12841
55	12790	11352	12838
50	12787	11351	12835
45	12784	11350	12832
40	12781	11349	12829
35	12778	11348	12826
30	12775	11347	12823
25	12772	11346	12820
20	12769	11345	12817
15	12766	11344	12814
15-40 Asst	12809	11359	13291
45-80 Asst	12810	11360	13292

## Reamers

Reamers are high quality, triangular bladed, stainless steel instruments designed to provide flexibility with integrity.



### Reamers (pkg of 6)

Size	21mm	25mm	30mm
45	06040	06132	
40	06038	06131	
35	06034	06127	
30	06033	06126	
25	06032	06120	
20	06031	06112	
15	06030	06111	
10	06029	06110	
08	14133	14134	
06	62252	62253	62254
15-40 Asst	821-2021	821-2025	