


Your practice is *our* inspiration.

Take 1 Advanced Designed for the way you work



Take 1 Advanced


A perfect fit for you, and your patients



Our expertise in high-performance impression materials, along with recent technological advancements, delivers a VPS material with the optimum combination of physical properties: strength, elasticity, dimensional stability, and the ability to register detail in any environment. Take 1 Advanced has the most extensive selection of set times, delivery systems and consistencies on the market today.

Research

As with any dental material, the physical and mechanical properties of impression material dictate how well it will perform in different clinical situations. While some materials focus on a single attribute, it's critical for these properties to work well together. Take 1[®] Advanced™ features individual properties that rate as good or better than any other VPS on the market today, with collective performance exceeding that of competitive products.



Features

Outstanding tear strength and elongation. Strong, flexible material stretches around undercuts without tearing or distorting.

Hydrophilic wash material. Displaces intraoral fluids to produce impressions with great detail and no voids or gaps.


Widest selection of viscosities, set times and delivery methods. Helps overcome common difficulties in impression taking, and matches doctor preference.



Internal testing. Data available upon request.

Tear Strength

Take 1 Advanced Wash material resists tearing around undercuts upon removal from the mouth. That combined with a high elongation percentage ensures preparation margins remain intact, even in tough cases below the gumline.




		LB Wash	RB Wash	Med/Mono	HB Tray	Rigid Tray	Putty	Bite Reg	
Options: Delivery, Vis and Set Times time	scosity	Flowable, highly flexible and hydrophilic	Highly thixotropic with incredible tear strength	The perfect balance between flow and body	Outstanding dimensional stability and elasticity	Designed to provide stability & support closed bite trays	Ultra-high viscosity for extra mechanical compression	Exceptional resistance to deformation to protect occlusal integrity	
Unidose® Kerr's unique Unidose de dual benefits: targeted s' delivery and the eliminat cross-contamination cor	yringe ion of	Super Fast Fast Set Regular Set	Super Fast Fast Set Regular Set						
Cartridge Convenient 50 ml cartrid mixing & application of Advanced precise and e	Take 1	Super Fast Fast Set Regular Set	Super Fast Fast Set Regular Set	Super Fast Fast Set Regular Set	Super Fast Fast Set Regular Set	Super Fast Fast Set		Fast Set	
Volume For busy offices, Take 1 / Volume delivery speeds of impression taking, reduce fatigue and lowers cost p	up es hand			Super Fast Fast Set Regular Set	Super Fast Fast Set Regular Set	Super Fast Fast Set			
Hand Mix For instances when a hi viscosity tray material is and other mixing metho are unavailable	desired						Regular Set		



Take 1 Advanced Bite Registration

Take 1 Advanced Bite completes the Take 1 Advanced portfolio providing you a complete solution within a single family of products.



Ordering

34064 Cartridge Intro Kit

Contains: 1 Cartridge (50 ml) LB Wash (Fast Set)
1 Cartridge (50 ml) RB Wash (Fast Set)
3 Cartridges (50 ml) Tray (Fast Set)
1 Cartridge (50 ml) Bite Registration

1 Extruder

6 Small Mix Tips

6 Intraoral Syringe Tips (Small Universal)

6 Large Mix Tips

Refills

Two Packs: 2 Cartridges (50 ml), 6 Mixing Tips

33953 Light Body Wash FS
33951 Light Body Wash RS
33955 Light Body Wash SF
33971 Regular Body Wash FS
33969 Regular Body Wash RS
33973 Regular Body Wash SF
(Wash viscosity 2-packs also contain intraoral tips)

 33959
 Medium FS

 33957
 Medium RS

 33961
 Medium SF

 33965
 Tray FS

 33963
 Tray RS

 33967
 Tray SF

 33967
 Tray SF

 34457
 Rigid Tray FS

 34458
 Rigid Tray SF

 33975
 Bite Registration

Value Packs: 24 Cartridges (50 ml)

33954 **Light Body Wash FS** 33952 **Light Body Wash RS** 33956 **Light Body Wash SF** 33972 Regular Body Wash FS Regular Body Wash RS 33970 Regular Body Wash SF 33974 33960 Medium FS 33958 Medium RS 33962 **Medium SF** 33966 Tray FS 33964 Tray RS 33968 Tray SF Rigid Tray FS 35186 35187 **Rigid Tray SF**

Unidose: 20 Syringes (2 ml)

33995 Light Body Wash FS 20 pack 33994 Light Body Wash RS 20 pack 33996 Light Body Wash SF 20 pack 33998 Regular Body Wash FS 20 pack 33997 Regular Body Wash RS 20 pack 33999 Regular Body Wash SF 20 pack

Volume: 2 Hard Cartridges (380 ml)

34695 Tray RS 34697 Tray SF 34696 Trav FS 34690 Medium/Monophase SF 34688 Medium/Monophase RS Medium/Monophase FS 34689 34825 Rigid Tray FS 34826 **Rigid Tray SF**

Putty: 1 Base (400 g), 1 Catalyst (400 g)

34070 Putty

Accessories

28747	Extruder
29363	Mix Tips, Small (pkg of 48)
29333	Mix Tips, Large (pkg of 48)
27517	Intraoral Syringe Tips, Small Univer
	(100 disposable)
20731	Intraoral Syringe Tips, Large
	(100 disposable)
16622	VPS Adhesive (15 ml bottle)
25777	VPS Adhesive (59 ml bottle)
35115	Volume Machine by Renfert
34746	Dynamic Mixing Tips for Volume
	(50 tips, 2 bayonet rings)

©2012 Kerr Corporation Lit. No. 34117D All trademarks are property of their respective owners.

Kerr