
KIMBERLY-CLARK[®] **PROTECTAPAD** Multi-purpose medical laminate

Kimberly-Clark PROTECTA**PAD** is made of 4 layers of paper tissue and 1 layer of soft polyethylene laminated together to provide a waterproof absorbent surface protection.

Used extensively in dental surgeries, Kimberly-Clark PROTECTA**PAD** is the ideal dental bib. It is used to provide a hygienic surface for dental tools and also for covering benches and trays in pathology labs.

Kimberly-Clark PROTECTA**PAD** is also used as a medical procedures underlay where fluids are involved, eg. Pillow protection.



PRODUCT ORDERING INFORMATION

Code	Product Description	Size	Sheets/ Pack	Case
2706	Kimberly-Clark PROTECTAPAD Small	28.5 x 21.5cm	8 packs of 10	800
2705	Kimberly-Clark PROTECTAPAD Large	28.5 x 43cm	4 packs of 100	400

	Kimberly-Clark
--	-----------------------

Kimberly-Clark Australia Pty Limited ABN 65 000 032 333 52 Alfred Street, Milsons Point NSW Australia 2061 To order: 1800 810 077 Health Care Advisory Service: 1800 101 021 Fax 1800 812 255 Kimberly-Clark New Zealand 70 Kerrs Road Manukau City Auckland New Zealand Phone: 09 262 8168

Fax: 09 262 8185 www.kca.com.au www.kchealthcare.com

Your Distributor:

® Registered Trademark Kimberly-Clark Worldwide Inc. © 2005 KCWW cfm/09.05/FA2

