



Ultra COMFORT, Ultra PROTECTION with an ERGONOMIC DESIGN

NEW



MICRO-TOUCH® Ultraform® UF-524: ERGONOMIC NITRILE EXAMINATION GLOVES

- Designed with **ERGOFORM™ Ergonomic Design Technology** to support musculoskeletal health
- Thin design for superior tactile sensitivity when handling instruments
- Soft nitrile formulation forms to hands for comfort during extended wear

ORDER INFORMATION

XS	[AN-UF-524-XS]	S	[AN-UF-524-S]	M	[AN-UF-524-M]
L	[AN-UF-524-L]	XL	[AN-UF-524-XL] (250 gloves / box)		

300 gloves / box to reduce waste

Now
Featuring



➤ You can order by contacting 0800 808 855 or learn more about the range of dentistry gloves by visiting henryschein.co.nz

MICR+FLEX®

Comfortable, thin nitrile examination gloves that offers affordable protection with an ergonomic design

MICROFLEX® Ultraform® UF-524

MICROFLEX® Ultraform® are thin nitrile examination gloves that reliably protect hands while offering exceptional tactility for instrument handling. Made with ERGOFORM™ Ergonomic Design Technology, these thin, tactile gloves are proven to reduce the amount of muscle effort dental practitioners exert while performing repetitive tasks when compared to similar disposable nitrile gloves¹.

PRODUCT INFORMATION

Material	Nitrile	Colour	Cobalt Blue
Cuff Length	Standard	External Glove Surface	Textured Fingertips
Palm Thickness mm (mils)	0.06 / 2.4	Watertight Barrier AQL	1.5
Finger Thickness mm (mils)	0.10 / 3.9	ISO Certification	ISO 13485:2012


ORDER INFORMATION

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)*
Henry Schein Product Code	AN-UF-524-XS	AN-UF-524-S	AN-UF-524-M	AN-UF-524-L	AN-UF-524-XL

PACKAGING AND STORAGE

Size	300 gloves per box; 10 boxes per case; 3000 gloves per case *XL size - 250 gloves per box; 2500 gloves per case
Shelf Life	3 years
Storage	Keep out of direct sunlight; store in a cool, dry place. Keep away from sources of ozone or ignition.

ERGOFORM™ Ergonomic Design Technology enables hand protection that supports musculoskeletal health during repetitive tasks and improves worker performance.



ERGOFORM™
Ergonomic Design Technology

1. Based on laboratory tests conducted by US Ergo and commissioned by Ansell.