

# DenJet®

SANITISER / CLEANER



**DS**

**WS**

**DS** - DAILY SANITISER / CLEANER  
**WS** - WEEKLY SANITISER / DESCALER

DenJet®-DS and DenJet®-WS are powerful effective foam-free clinical detergent/sanitiser. Formulated to clean & sanitise dental aspiration/suction lines & separation systems.

## **DS** DenJet®-DS (Daily)

Use once or twice per day as part of the Practice Management instructions for the maintenance & sanitation of aspiration lines for assisting the removal of bioburden.

An effective sanitiser is used that incorporates an advanced de-foaming formulation for rapid sanitation of the aspiration lines. The formula incorporates unique de-foaming surfactants to ensure detergency of contaminating soils and the penetration of the biofilm. Mildly alkaline auxiliary agents enhance the activity of the sanitising agents and will aid in the solubilization of proteins and soils.

## **WS** DenJet®-WS (Weekly)

DenJet®-WS should be used once or twice per week as part of the routine maintenance of aspiration lines to assist the removal of scale and biofilm.

A Highly effective descaler that will help to remove soils that are resistant to other detergents. Incorporates special de-foaming surfactants that will ensure the penetration of the biofilm and contaminating soils.

Mild acid descaling agents based on natural lactic acid and citric acid are incorporated for their excellent compatibility with metal, plastic and rubber components.

**Powerful effective  
foam-free aspiration  
& suction line cleaner**

### PRODUCT CODES:

DenJet®-DS 5 Litre Container - DJDS5  
DenJet®-DS 15 Litre Container - DJDS15  
DenJet®-WS 1 Litre Container - DJWS1



Specially formulated to be suitable for use with dental aspiration/suction pumps and lines.

Also available from:



Marketed by:  
**Majac Medical Products Pty. Ltd.**  
Tel: 1300 138 578 • Int: +61 7 3265 6355  
Email: sales@majacmedical.com.au  
[www.majacmedical.com.au](http://www.majacmedical.com.au)



FS 540017

# DenJet<sup>®</sup>

SANITISER / CLEANER

## Preventative Maintenance of **DENTAL ASPIRATION SYSTEMS**

**DS** **WS** **DS** - DAILY SANITISER / CLEANER  
**WS** - WEEKLY SANITISER / DESCALER

Routine cleaning is an essential part in keeping the aspiration lines in working condition and a guarantee of effective trouble free operation.

Sanitation of the waste system is essential to maintaining a hygienic aspiration system that will protect patients and staff from hazardous microorganisms.

**IMPORTANT:**

Strong domestic hypochlorite solutions should never be used to clean aspiration lines because the chlorine can damage rubber and metal components. These types of cleaners are usually caustic and therefore can damage some materials, they are also high foaming.

Phenolic compounds in some industrial disinfectants will damage plastics and rubber.

**DenJet<sup>®</sup>-DS is a use daily mild alkaline sanitiser/cleaner.**

DenJet<sup>®</sup>-DS will remove protein and bio-soils from the suction tubes while eliminating bacteria before they can form bio-films that will narrow the lumen of the lines and straining the pumps.

**A weekly descale with DenJet<sup>®</sup>-WS will remove material that is resistant to other cleaners.**

A 'NO' foaming sanitiser/cleaner is important to maintain the efficiency of the pumps and traps. Cleaning and Sanitising on a daily basis prevents the build-up of bacterial growth (biofilm) which can trap other debris.

