

MAXIMA® SELF SEAL
**STERILISATION
POUCHES**

SECURELY SEALED.
**STERILISATION
MAINTAINED.**



Order today at [HENRYSCHEIN.CO.NZ](https://henryschein.co.nz)



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS

MAXIMA® SELF SEAL STERILISATION POUCHES

- Dual indicators change color when a sterilisation process has begun. To confirm sterilisation, the use of a biological monitor is necessary
- Steam Autoclaves: Turns from blue to black at 121°C for a minimum of 20 minutes
- Ethylene Oxide (EtO): Turns from pink to golden brown after 60 minutes at 50°C, RH 50–80%, and EO concentration of 735 mg/L
- Reliably meets sterilisation standards when used according to the Instructions for Use (IFU)
- When properly stored, instruments maintain their sterility for up to 1 year post sterilisation
- Compatible with autoclaves and ethylene oxide (EtO) chamber sterilisers



STEAM STERILISATION



EO STERILISATION



SELF SEAL STERILISATION POUCHES 200/pack

ITEM CODE	SIZE (MM)	SIZE (INCHES)
(9006153)	57mm x 100mm	2.25x4"
(9006154)	70mm x 230mm	2.75x9"
(9006155)	90mm x 140mm	3.5x 5.5"
(9006156)	90mm x 255mm	3.5x10"
(9006157)	133mm x 255mm	5.25x10"
(9006158)	190mm x 330mm	7.5x13"
(9006159)	230mm x 380mm	9x15"