



Opalescence™ Whitening Storage Guide

	Room Temperature Storage	Refrigeration Required	Shelf Life
Opalescence Go™ Whitening Trays Prefilled		✓	24 months
Opalescence™ Boost™ 40% Office Whitening In-	✓		24 months
Opalescence™ Boost™ 35% Office Whitening In-		✓	24 months
Opalescence™ 10% PF Take-Home Whitening Gel Non-		✓	24 months
Opalescence™ PF 10%, 15%, 20% Take-Home Whitening Gel	✓		24 months
Opalescence™ PF 35%, 45% Take-Home Whitening Gel		✓	24 months
Opalescence™ Endo Whitening Gel		✓	24 months
Opalustre™ Microabrasion Slurry	✓		24 months
UltraEZ™ Gel Syringes Desensitizing	✓		48 months
UltraEZ™ Gel UltraFit™ Trays Desensitizing	✓		18 months
OpalDam™ and OpalDam Green™ Light-Cured Resin Barriers	✓		48 months
Opalescence Office™ Whitening		✓	24 months

1800 29 09 29 | info.anz@ultradent.com

Ultradent Products Australia Pty Ltd PO Box 184 Manly, NSW 1655 Australia